



## Molecular cloning and genetic engineering – Competent Cell Production

### ● Aim

Competent cell: physical and chemical methods induce the cells to be in the optimum physiological state of uptake and accommodation of foreign DNA.

### ● Materials

0.1mol/L  $\text{CaCl}_2$

LB plate

ddH<sub>2</sub>O

### ● Procedure

1. Ethanol and ultraviolet treat all working areas for sterility.
2. Streak DH5 $\alpha$  cells on an LB plate and grow for single colonies at 37°C;
3. Fill an ice bucket halfway with ice. Use the ice to pre-chill centrifuge tubes and 0.1mol/L  $\text{CaCl}_2$ ;
4. Pick a single colony in 30 ml LB medium, shake it at 37°C, 160 r/min, until  $\text{OD}_{600} = 0.3$ , then take it out and put it on ice for 15minutes
5. Pipette 1ml of bacterial solution and sterilized in 1.5 ml centrifugal tube.



---

Centrifugation at 4°C for 5min at 5000r/min;

6. Discarded the supernatant, suck up the residual culture medium and gather cells.

7. Add 500 ul pre-chilled 0.1mol/L  $\text{CaCl}_2$  to resuspend the cells, ice bath for 20 min.

8. Centrifugation at 4°C for 5 min at 5000r/min;

9. Discarded the supernatant, suck up the residual water. Add 100 ul pre-chilled 0.1mol/L  $\text{CaCl}_2$ , place at 4°C for confirmation.

