

DNS Glucose Assay

Introduction

The procedure determine the glucose concentration of a sample using 3,5-dinitrosalicylic acid (DNS). Adapted by Philip Sørensen and Jacob Mejlsted from [D. Navarro et al.](#) and [Worthington biochem.](#)

Materials

› Materials

- › Plate reader

› DNS reagent

- › 1 grams of 3,5-dinitrosalicylic acid
- › Milli-Q water
- › 30.0 grams sodium potassium tartrate tetrahydrate
- › 20 mL 2 M NaOH

› Consumables

- › Microtiter plate

Procedure

Preparation of DNS reagent

1. Dissolve 1 g of 3,5-dinitrosalicylic acid in 50 mL Milli-Q water
2. Add slowly 30.0 grams of sodium potassium tartrate tetrahydrate
3. Add 20 mL of 2 M NaOH
4. Dilute to a final volume of 100 mL with reagent grade water. Protect from carbon dioxide and store no longer than 2 weeks.

Running the assay

5. 100 μ L media/sample is added to the microtiter plate
6. 100 μ L DNS reagent is added to each sample
7. Plate is heated to 98 $^{\circ}$ C for 10 min
8. Absorbance at 540 nm is measured
9. **Remember to include a standard to determine actual concentrations**