



## COLOR CODING KEY

**Characterization**

**Construct**

**Data Collection/Analysis**

**Improvement of Biobrick**

**Plant-Care**

**Plasmid**

**Plasmid/Construct Design**

**Preparatory Work**

---

### **Tuesday April 30th**

- Received 300 seeds from Martin Lab

---

### **Friday May 3rd**

- Planted first batch of seeds as practice

---

### **Wednesday May 29th**

- Two new batches of seeds planted

---

### **Friday June 7th**

- Received seeds from Tran Lab
- Planted four new batches of seeds

---

### **Thursday June 20th**

- Started planning for wet lab to start tomorrow
- LB powder, autoclave (find a working autoclave), make sterile LB media, put in autoclave/cool (sterilize media), aliquot into tubes, stab bacteria and put in liquid media (not airtight), put in shaker

#### *Making LB media*

- 25 g of powder for every 1 L of water
- Autoclave for 15 min on 120 degrees Celsius
- Prepare hot water bath for effective cooling

- Allow to cool before making additions (including antibiotics)
- 

### Friday June 21st

- All wet lab materials came in
- Constructs and primers are in drawer at the front
- Backbones, knockout yeast, and Rowley plasmid are in 4 degree fridge

### Made agar plates

- 25g of LB powder, 12 g of agar powder
  - Add water, and autoclave
  - Sterilize bench/cool in water bath
  - Add antibiotics (Ampicillin)
-