

Name: Kennex, Laura, Asma, Amirah, Sijia, Jaizi, Krithika, Chiara

Date: 7/19/19

Goals:

1. Overnight cultures
 - a. Pcb302 in A. Tume from transformations done on 7/17/19
2. Restriction Digest
 - a. DinolIII-P2 digested with BglIII and XbaI
3. Gel electrophoresis
 - a. DinolIII P2 digested with BglIII & XbaI
 - b. Codon optimized RFP digested with BglIII & XbaI

Name: Chiara

Date: 7/19/19

Goal:

1. Overnight cultures
 - a. Pcb302 in A. Tumefaciens from transformations done on 7/17/19

Protocol:

Overnight Cultures

1. Added about 7 mL of LB with 7 μ L kanamycin to a 15 mL Falcon tube
2. Dipped a p10 tip into the selected colonies and drop into the tube
3. Incubated in the water bath at 28° C at 200 rpm over the weekend

Results:

N/A

Conclusion:

N/A

Date: 07/19/2019

Goals:

1. Restriction Digest
 - a. DinIII-P2 samples digested with BglII and XbaI

Protocol:


30 μ L Fast Digest Restriction Digest

1. Prepared a Fast Digest concentration cocktail with the following proportions: 1 μ L XbaI, 1 μ L BglII, 3 μ L of 10X Fast Digest Buffer, and 15 μ L of diH₂O.
2. Added 20 μ L of this cocktail to a clean 1.5 Eppendorf tube and then added 10 μ L of DNA
3. Incubated at 37° C for 30 minutes.

Results:

[illegible]

Expected Gel-



MW: GeneRuler™ 1 kb Plus DNA Ladder

2:

3:

4:

5:

6:

7:

8: DinoIII Part 2 in pUCIDT vector

BglII + XbaI

1. 4612 bp = BglII (2240) - XbaI (1514)

2. 726 bp = XbaI (1514) - BglII (2240)


9: DinoIII Part 2 in pUCIDT vector

XbaI + BglII

1. 4612 bp = BglII (2240) - XbaI (1514)

2. 726 bp = XbaI (1514) - BglII (2240)

Observed gel-



Conclusion:

The digest was successful; the observed bands are similar to the expected bands. The next step is gel extraction.

Name: Amirah

Date: 7/19/19

Goal:

1. Ran RFP codon optimized pre-done restriction digest on gel

Protocol:

Preparing, Loading, and Running a 1% Agarose Gel

Preparing

1. Added 1 g of Agarose in 100 mL of 1X TBE in an Erlenmeyer flask
2. Heated in the microwave until fully dissolved
3. Allowed the solution to cool until comfortable to touch
4. Added 10 μ L GelRed Nucleic Acid Gel Stain and mixed
5. Inserted casting tray, made sure the rubber on the sides was not overlapping
6. Carefully poured the agarose into the tray and placed the comb to create the wells
7. Allowed the gel to solidify
8. Once solidified, changed the orientation of casting tray where the rubber sides are not in contact with the sides of the system.
9. Poured in 1X TBE into the gel electrophoresis system to the fill line, being sure to submerge the gel, and removed the comb

Loading

1. Loaded ~5 μ L of the ladder in the first well
2. Prepared your samples to load by adding in 1 μ L of 6X Loading dye for every 5 μ L of DNA and loaded
 - a. 10 μ L of DNA was loaded

Running

1. Once the gel had been loaded, slid on the cover making sure the negative electrode is closest to the DNA and the positive electrode is at the bottom of the gel
2. Ran for about an hour at 150 V

Results:

Gel did not look like simulation

Conclusion:

Did not do extraction. Left gel in fridge