

Name: Kennex, Rehmat, Asma, Yujie Huang, Jiayi Lan, Xuecheng Ye

Date: 8/2/19

Goal:

1. Run a gel on the PCR products from yesterday
2. Run a gel and do gel extraction on the digested codon optimized RFP
3. Ligation of RFP and DinIII plasmid
4. Ligation for composite part

Name: Yujie Huang, Jiayi Lan, Xuecheng Ye

Date: 8/2/19

Goal:

1. Run a gel on the PCR products from yesterday

Protocol:

## **Preparing, Loading, and Running a 1% Agarose Gel**

### **Preparing**

1. Added 1 g of Agarose in 100 mL of 1X TBE in an Erlenmeyer flask.
2. Heated until fully dissolved.
3. Added 10  $\mu$ L SYBR Gold Nucleic Acid Gel Stain when it cooled enough to touch.
4. Inserted casting tray.
5. Poured the agarose into the tray and placed the comb to create the wells. Gel solidified
6. Changed the orientation of casting tray so the rubber sides were not in contact with the sides of the system.
7. Poured in 1X TBE into the gel electrophoresis system to the fill line, making sure to submerge the gel.

### **Loading**

1. Loaded 5  $\mu$ L of the GeneRuler 1kb Plus ladder in the first well .
2. Prepared samples to load by adding in 1  $\mu$ L of 6X Loading dye for every 5  $\mu$ L of DNA and loaded

### **Running**

1. Ran for 1 hour at 100 volts.

Results:




Figure 1: from right to left is ladder, 1-12




Figure 2: from left to right is ladder, 13-24

Name: Kennex, Asma

Date: 8/2/19

Goal:

1. Gel Extraction of RFP Midiprep #3 and #7 from 7/24 and RFP Miniprep #10 from 7/11

Protocol:

## **Preparing, Loading, and Running a 1% Agarose Gel**

### **Preparing**

1. Added 1 g of Agarose in 100 mL of 1X TBE in an Erlenmeyer flask.
2. Heated until fully dissolved.
3. Added 10  $\mu$ L GelRed Nucleic Acid Gel Stain when it cooled enough to touch.
4. Inserted casting tray.
5. Poured the agarose into the tray and placed the comb to create the wells. Gel solidified
6. Changed the orientation of casting tray so the rubber sides were not in contact with the sides of the system.
7. Poured in 1X TBE into the gel electrophoresis system to the fill line, making sure to submerge the gel.


### **Loading**

1. Loaded 10  $\mu$ L of the GeneRuler 1kb Plus ladder in the first well .
2. Loaded 30  $\mu$ L of digested RFP #3 + 6 uL of Purple Loading Dye into third well.
3. Loaded 30 uL of digested RFP #7 + 6 uL of Purple Loading Dye into fifth well.
4. Loaded 30 uL of digested RFP #10+ 6 uL of Purple Loading Dye into seventh well.

### **Running**

1. Ran for 1 hour at 100 volts.

## Results:



### **Gel Extraction**

QIAQuick Gel Extraction Kit Lot # 42141174

1. Ran a restriction digest on the targeted DNA part using restriction enzymes yesterday and ran an agarose gel for 1 hour before starting
2. Cut the targeted DNA sequence out using a razor blade, making sure to get as much DNA while limiting the amount of agarose extracted
3. Pre-weighed empty Eppendorf tubes before adding the gel excisions.
4. Added the gel extracts to Eppendorf tubes and weighed again.
5. Calculated the mass of the gel using the difference of the two measurements.
6. The RFP fragment 3 weighed 350 mg, RFP fragment 7 weighed 330 mg, and the RFP fragment 10 weighed 400 mg.
7. Added 1050  $\mu$ L of Buffer QG to the RFP fragment 3, 990  $\mu$ L to the RFP fragment 7, and 1200  $\mu$ L to the RFP fragment 10.
8. Incubated the tubes at 50° C for 10 minutes and vortexed every 2 minutes to help dissolve the gel
  - a. Checked to make sure the color of the mixture is yellow
9. Once dissolved, added 350  $\mu$ L of isopropanol to the RFP fragment 3, 330  $\mu$ L to the RFP fragment 7, and 400  $\mu$ L to the RFP fragment 10, and mixed.
10. Placed a QIAquick spin column in a provided 2 mL collection tube.
11. Added 700  $\mu$ L of the solution to the spin column at a time and centrifuged at 13,000 rpm for 1 minute and discarded the flow through. Repeated until all of the solution had ran through.
12. Added 500  $\mu$ L of Buffer QG to the spin columns to remove traces of agarose and centrifuged for 1 minute.
13. Added 750  $\mu$ L of Buffer PE to the column to wash and centrifuged for 1 minute.

14. Discarded the flow through and centrifuged for an additional 1 minute at 13,000 rpm to remove residual buffer.
15. Placed the spin column in a clean 1.5 mL Eppendorf tube
16. Added 40  $\mu\text{L}$  of warmed Buffer EB to the center of the spin column, allowed to sit for 1 minute, and centrifuged for 1 minute.
17. Measured and recorded the concentrations.

### Results

RFP 3	12.5 ng/ $\mu\text{L}$
RFP 7	10 ng/ $\mu\text{L}$
RFP	5 ng/ $\mu\text{L}$

Lane 1: 1 KB Plus Ladder

Lane 2: Empty


Lane 3: RFP 3

Lane 4: Empty

Lane 5: RFP 7

Lane 6: Empty

Lane 7: RFP 10



Name: Rehmat, Kennex, Asma

Date: 8/2/19

Goal:

1. Ligate the DinIII plasmid and the codon optimized RFP

Protocol:

Ligation

1. Added 6  $\mu\text{L}$  of  $\text{dH}_2\text{O}$  to a clean 1.5 mL Eppendorf tube, 1  $\mu\text{L}$  of T4 DNA Ligase Buffer, 1  $\mu\text{L}$  of the Dino III plasmid and 1  $\mu\text{L}$  of the codon optimized RFP part and mixed.
2. Added 1  $\mu\text{L}$  of T4 DNA Ligase.
3. Pipet mixed the tube and incubated at room temperature for 10 minutes.

Conclusion:

We will transform this ligation into DH5alpha E. Coli cells to see if our ligation worked.