

Name: Laura, Beenle

Date: 10/17/19

Goal: Transform Symbiodinium and O. Marina using Lonza

Protocol:

1. 1.5 mL of the symbiodinium (cell count: 1,092,000 cells/mL) was added to an eppendorf tube and spun down at 900g for 3 minutes. 3.5 mL of the o. Marina (cell count: 420,000 cells/mL) was added to a 50mL falcon tube and spun down at 1000g for 3 minutes.
2. The supernatant was removed from both tubes and 98 μ L of SG and 22 μ L of the supplement were added to EACH of the tubes and were swirled to resuspend the pellet.
3. 40 μ L aliquots were made of each cell type (6 tubes total)
4. O. marina tubes
 - a. 20 μ L of Dino III RFP (ng/ μ L)
 - b. 20 μ L of Dino III GFP (ng/ μ L)
 - c. 20 μ L diH₂O (for the blank)
5. S. microadriaticum tubes
 - a. 40 μ L of Dino III RFP (ng/ μ L)
 - b. 20 μ L of Dino III GFP (ng/ μ L)
 - c. 20 μ L diH₂O (for the blank)
6. The Lonza machine was set for the correct pulse codes (DS-137 and DS-130) and 25 μ L of each sample was added to the corresponding cuvette well
7. After shocking, 80 μ L of media (ASP8A for s. microadriaticum and F/2 for o. marina) was added to the cuvette well and pipetted up and down three times. The contents of the cuvette were then transferred to the corresponding well in a 24-well plate.
8. The 24-well plates were placed in artificial light (hub) for incubation.

Notes:

- On the 24-well plate, A stands for pulse code DS-137 and B stands for pulse code DS-130
- The 24-well plate on top contains transformed O. marina and the one on the bottom contains transformed S. microadriaticum