

# USER cloning

---

## Introduction

This is a protocol to do USER cloning based on the standard PacI/Nt.BbvCI cassette plasmid

It's adapted from protocols by Rasmus Frandsen and Kresten J. Olsen (DTU Bioengineering) respectively

## Materials

### › Reagents

- › USER-linearized insert(s)
- › Linearised plasmid backbone (Plasmid digested with PacI/Nt.BbvCI)
- › USER enzyme
- › Cutsmart buffer
- › Milli-Q water

### › Consumables

- › PCR strips

## Procedure

### Preparing USER-reaction

1. In a PCR-tube, combine;
  - 200 ng USER-linearized insert(s),
  - 100 ng USER-linearized vector backbone,
  - 1 uL 10X CutSmart buffer,
  - 1 uL USER enzyme,
  - H2O to 10 uL
2. Incubate reaction mix as described below
  - 37C for 25 min
  - 25C for 10 min
  - 20C for 10 min,
  - 15C for 10 min
  - 10C until further use
3. Transform reaction mixture in E. coli as described in the Hifi assembly protocol