

Name: Chiara Brust

Date: 9/19/19

Goal:

1. Transform *S. Microadriaticum* using fungi electroporation protocol

Name: Chiara Brust

Date: 9/19/19

Goal:

1. Transform *S. Microadriaticum* using fungi electroporation protocol

Protocol:

1. Cultured *Symbiodinium Microadriaticum* cells in ASP-8A with filtered seawater medium under natural sunlight for about 2 months and *O. Marina* cells in f/2 media for about 2 months with no antibiotics

- Cultures:
  - *S. Microadriaticum*: ASP 8a-SW media 100 mL labelled 8/12/19; [cell]=  $1.0 \times 10^5$  cells/mL

2. Harvested the cells by centrifugation at 800 g for 5 min at 4°C.

- a. *S. Microadriaticum*: Pelleted 10 mL 3 samples =  $1.0 \times 10^6$  cells/mL

3. Used 500  $\mu$ L of 0.1M EDTA to resuspend the cell pellet

- a. [S. Micro]=  $2.58 \times 10^6$  cells/mL
- b. [S. Micro]=  $2.78 \times 10^6$  cells/mL
- c. [S. Micro]=  $2.00 \times 10^6$  cells/mL

4. Centrifuged at 800 g for 2 min at 4°C.

5. Washed cells with 500  $\mu$ L of 10% Glycerol 3 times, centrifuged at 800 g for 2 min in 4°C.

6. Resuspended the pellet in 50  $\mu\text{L}$  of 10% Glycerol
  - a. [S. Micro]=  $1.18 \times 10^6$  cells/mL
  - b. [S. Micro]=  $5.2 \times 10^5$  cells/mL
  - c. [S. Micro]=  $1.32 \times 10^6$  cells/mL
8. Incubated  $\sim 40\mu\text{L}$  of cells with no DNA on ice for 5 min.
9. Put cells into a 0.2 cm cuvette, and electroporated using SHS (2.0 kV, 1 pulse), SC2 (1.5 kV, 1 pulse), or DIC (1.0 kV, 2 pulses, 1.0 msec) program with Bio-Rad MicroPulser 165-2100.
10. Added 1mL of ASP-8A SW medium to the 0.2 cm cuvette, mixed well, and transferred to a 15 mL tube
  - a. [S. Micro]=  $8.0 \times 10^4$  cells/mL
  - b. [S. Micro]=  $1.2 \times 10^5$  cells/mL
  - c. [S. Micro]=  $4.0 \times 10^4$  cells/mL
11. Left the cultures to grow under natural sunlight
12. Observed the cells under a normal microscope in 1-3 days and according to the need.