

Name: Sijia, Jaizi, Rehamt, Chiara, Shakera, Kennex, Amirah

Date: 7/15/19

Goals:

1. Make glycerol stock for overnight cultures for the transformations done on 7/9/19
 - a. Pcb302 in A. Tumefaciens from papers 1 & 2
2. Miniprep for pcb302 in A. Tume from papers 1 & 2 overnight cultures for the transformations done on 7/9/19
3. Transform
 - a. K1357009
 - b. Dinolll part 2
4. Run gel of K1357009 colony PCR
5. Overnight cultures
 - a. PGEX-HCG from glycerol stocks

Name: Rehmat Babar, Shakera Thomas

Date: 7/15/19

Goal:

1. Miniprep Extract pCB302-gfp-MBD plasmid from agrobacterium overnights

Protocol:

Mini Preps for *Agrobacterium tumefaciens*

1. Centrifuged 10 mL of overnight for 15 minutes at 3500 rpm and resuspended in 250 μ L buffer P1 containing 0.1 mg/mL RNase A.
2. Added 250 μ L lysis buffer P2 to the tube and inverted gently 6 times to mix.
3. Added 350 μ L neutralization buffer N3 to the tube and inverted immediately but gently 6 times.
4. Centrifuged the lysate for 10 min at 13,000 rpm
5. Placed a QIAprep Spin Column in a 2 mL collection tube.
6. Transferred the cleared lysates from step 4 to the QIAprep Spin Column by decanting or pipetting.
7. Centrifuged for 60 seconds at 13,000 rpm and discarded flow-through.
8. Washed the QIAprep Spin Column by adding 500 μ L of Buffer PB and centrifuged for 60 seconds at 13,000 rpm and discarded flow-through.
9. Washed the QIAprep Spin Column by adding 750 μ L of Buffer PE and centrifuging at 60 seconds at 13,000 rpm and discarded the flow through.
10. Centrifuged for an additional 1 min to remove residual wash buffer at 13,000 rpm.
11. Placed the QIAprep Spin Column in a clean 1.5 mL microcentrifuge tube.
12. Added 50 μ L of Buffer EB to the center of each QIAprep Spin Column, let stand for 1 min, and centrifuged for 1 min.

Results:

Samples	DNA Concentrations (ng/ μ L)
1	0.35
2	0.15
3	0.40
4	0.25
5	0.20
6	0.10

7	0.20
8	0.40
9	0.10
10	0.25
11	0.15
12	0.15
13	0.25
14	0.40
15	0.40
16	0.40
17	0.10
18	0.10
19	0.15
20	0.15
21	0.15
22	0.20
23	0.15
24	0.25
25	2.00
26	0.10
27	0.10
28	0.15
29	0.15
30	0.10
31	0.25
32	0.10

33	0.15
34	0.55
35	0.55
36	0.15
37A*	0.25
37B*	0.55
38	0.60
40	0.50
41	0.35
42	N/A
43	0.15
44	0.10
45	0.15
46	0.20
47	0.15
48	0.20
49	0.20
50	0.15
51	0.90
52	0.35
53	0.45

Conclusion:

The DNA concentrations were too low

*These samples were mislabeled and consequently unidentifiable.

Name: Sijia Qin, Jiazi Tian

Date: 7/15/19

Goals:

1. Make glycerol stock for overnight cultures for the transformations done on 7/9/19
 - a. Pcb302

Protocol:

1. 500 μ l of 53 samples and 500 μ l of 50% glycerol was added.
2. The box was labeled as glycerol stock of o/n cultures on 7/9/19 in the top of freezer.

Conclusion:

N/A

Name: Chiara Brust

Date: 7/15/19

Goal:

1. Transform K1357009

Protocol:

Heat Shock Transformation

1. Thawed One Shot TOP10 chemically competent cells on ice.
2. Added 2 μL of DNA sample into competent cells
 - a. K1357009
 - b. RFP positive control
3. Incubated the cells on ice for 35 minutes.
4. After the ice incubation, placed the samples into a 42° C water bath for 30 seconds.
5. **Quickly** took them out and **immediately** added 250 μL of SOC medium
6. Placed the samples into a 37° C shaking water incubator for 1 hour at 200 rpm.
7. After shaking for 1 hour, smeared 150 μL of the solution onto an agar plate with the respective antibiotics.
 - a. Chloramphenicol
8. Incubated plates at 37°C for at least 24 hours.

Results:

N/A

Conclusion:

N/A

Name: Chiara Brust

Date: 7/15/19

Goal:

1. Resuspend DinolIII Part 2 plasmid
2. Transform DinolIII part 2 plasmid into E. Coli

Protocol:

Resuspension

1. Centrifuged the tube prior to opening at 13,000 RPM for 5 minutes
2. Resuspended DNA in IDTE buffer to the desired concentration
3. Vortexed for 20 seconds
4. Incubated the tube at room temperature for 30 minutes
5. Centrifuged for 1 minute

Transformation

Heat Shock

1. Thawed One Shot TOP10 chemically competent cells on ice.
2. Added 1 μL of DNA sample into competent cells
 - a. DinolIII part 2
3. Incubated the cells on ice for 35 minutes.
4. After the ice incubation, placed the samples into a 42° C water bath for 30 seconds.
5. **Quickly** took them out and **immediately** added 250 μL of SOC medium
6. Placed the samples into a 37° C shaking water incubator for 1 hour at 200 rpm.
7. After shaking for 1 hour, smeared 150 μL of the solution onto an agar plate with the respective antibiotics.
 - a. Ampicillin
8. Incubated plates at 37°C for at least 24 hours.

Results:

N/A

Conclusion:

N/A

Name: Chiara Brust

Date: 7/15/19

Goal:

1. Overnight cultures of PGEX-HCG from glycerol stocks

Protocol:

Overnight Cultures

1. Added about 6 mL of LB to a 15 mL Falcon tube along with 6 μ L of antibiotics
 - a. 1000X Chloramphenicol
2. Scraped some of the glycerol stock ice with the p10 tip and dropped into the tube
3. Incubated at 37° C at 220 rpm for 16-18 hours

Results:

N/A

Conclusion:

N/A

Name: Amirah Hurst

Date: 7/15/19

Goals:


1. To run gel electrophoresis on the PCR product from the K1357009 plasmid done on 7/12/19
 - a. Samples:
 - i. A1
 - ii. A2
 - iii. A3
 - iv. B1
 - v. B2
 - vi. B3
 - vii. C1
 - viii. C2
 - ix. C3
 - x. D1
 - xi. D2
 - xii. D3

Protocol:

1. Made two 1% Agarose gel
 - a. 1 g Agarose
 - b. 100 ml 1x TBE
 - i. Mixed in Erlenmeyer flask
 - ii. Heat in microwave until fully dissolved, swirling flask intermittently
 - c. Solution cooled
 - d. Added 10 ul of gel red
 - e. Made gel with cast
2. Loaded 5 ul of samples and ladder onto gel
3. Ran gel at 130v at high voltage for one hour


Results:

Gel 1



- Gel 1 key
 1. 1 kb ladder
 2. 100 bp ladder
 3. A1
 - 4.
 5. A2
 - 6.
 7. A3
 8. B1
 9. B2
 10. B3

Gel 2



- Gel 2 key
 1. 1 kb ladder
 - 2.
 3. C1
 - 4.
 5. C2
 - 6.
 7. C3
 - 8.
 9. D1
 - 10.
 11. D2
 - 12.
 13. D3
 14. 100 bp ladder