

Name: Asma, Krithika

Date: 9/25/19

Goal:

1. Transform pcb302 into *A. Tumefaciens*
 - a. Use DNA extracted from *E. Coli* transformations Plate B Colony 4
2. Test (blank) yeast transformation protocol on *S. microadriaticum*

Name: Asma

Date: 9/25/19

Goal:

3. Transform pcb302 into A. Tumefaciens
 - a. Use DNA extracted from E. Coli transformations Plate B Colony 4

Agrobacterium tumefaciens LBA4404

Protocol:

Electroporation of Agrobacterium tumefaciens

1. Thawed Agrobacterium tumefaciens cells on wet ice
2. Combined 1 μ L of pCB302-gfp-MBD plasmid DNA and 20 μ L of cells in an Eppendorf Tube
3. Pipetted the cells into a cuvette and electroporated at 2 kV
4. Added 1 mL of YM media and transferred to a 15 mL falcon tube
5. The tubes were incubated at 30°C at 200 rpm for 3 hours
6. 400 μ L of each culture was streaked onto a LB kanamycin plate.
7. 300 μ L of cultures 1 & 3 was streaked onto a YM kanamycin plate.
8. 300 μ L of culture 2 was streaked onto a LB kanamycin plate.
9. 200 μ L of each culture was also streaked onto a LB kanamycin plate.
10. The plates were incubated at 30°C for 48 hours

*Incubating 9/25 @6:30pm

Name: Krithika

Date: 9/25/19

Goal:

4. Blank transformation on *S. microadriaticum* using yeast protocol

Protocol:

1. 90 μ L of symbiodinium at an OD of 0.678 was added to an ice cold electroporation cuvette
2. The symbiodinium was incubated on ice for 5 minutes
3. The cells were pulsed at 1.5V and given 1mL 10% glycerol recovery
4. The cells were transferred to a 15 mL tube and placed in a non-shaking water bath at 30 °C for 1 hour
5. The cells were then put in 2 mL of ASP8A media and placed in artificial light for rapid growth