

Name: Saleh

Date: 8/23/19

Goal:

1. Gel extraction of Dino III RFP

Name: Saleh

Date: 8/23/19

Goal:

2. Gel extraction of Dino III RFP

Protocol:

Gel Extraction Dino III RFP

1. Used 50 mL falcon tube for each gel and added all gel fragments from one gel into the same 50 mL falcon tube
2. Weighed the gel fragments
3. Added 3 times volume of Buffer QG to gel fragments
 - a. Ex: Dissolved 3.6 grams of agarose in 11 mL Buffer QG
4. Vortexed to dissolve gel fragments
5. Added 1 volume of isopropanol
 - a. 3.6 mL isopropanol
6. Added 750 μ L of the solution to spin columns, centrifuged for 1 minute and discarded the flow through until all the solution had run through
 - a. For every 0.4 grams of agarose you need 1 spin column so for 3.6 g of agarose you'd need 9 spin columns
7. Added 750 μ L of Buffer QG to each of the spin columns, centrifuged for 1 minute, and discarded the flow through
8. Added 750 μ L Buffer PE to each of the spin columns, centrifuged for 1 minute, and discarded the flow through
9. Centrifuged the spin columns empty to remove the residual buffer.
10. Added 50 μ L Buffer EB to each of the spin columns and eluted the DNA into two different clean eppendorf tubes, making sure the solution doesn't touch the tip of the spin column.

Results:

The concentrations were too low to read.

Conclusion:

Since the concentrations were too low, an ethanol precipitation was subsequently prepared.