

Name: Asma, Kennex, Krithika, Chiara

Date: 08/07/2019

Goals:

1. Transform mCherry into BL21 and NEBa cells
 - a. Repeat transformation of mCherry with control
2. PCR on DinIII-rfp colonies from transformations done on 8/6/19
3. Restriction Digest codon optimized RFP 200 μ L reaction
 - a. Enzymes: XbaI & BglII
4. Overnight Culture on Ligation of DinIII and RFP

Name: Chiara Brust

Date: 8/7/19

Goal:

1. Restriction Digest codon optimized RFP 200 μL reaction
 - a. Enzymes: XbaI & BglII
 - b. Used C.O. RFP sample # 3 from 7/24/19 midiprep

Protocol:

Restriction Digest

60 μL Fast Digest Restriction Digest

1. Prepared a Fast Digest concentration cocktail with the following proportions: 2 μL Restriction Enzyme XbaI, 2 μL Restriction Enzyme BglII, 6 μL of 10X Fast Digest Buffer, and 30 μL of dH_2O .
2. Add 40 μL of this cocktail to a clean 1.5 Eppendorf tube and then add 20 μL DNA
3. Incubate at 37° C for 1 hour

Name: Kennex Lam and Krithika Karunakaran

Date: 8/7/19

Goal:

1. Overnight on Ligation of DinIII and RFP

Protocol:

1. 15 mL of LB + 15 uL of Ampicillin was added into 6 falcon tubes.
2. An individual colony was picked and a tip was dropped in a separate tube.
3. The tubes were left to shake at 220 rpm at 37° C overnight.

Results:

Conclusion:

We are going to measure to concentration of DinIII with RFP and perform a restriction digest on the plasmid in order to linearize it before transforming the ligated DinIII into *O. marina*.

Name: Asma Khimani

Date: 08/07/2019

Goals:

1. PCR Colony on *DinIII-rfp* colonies from transformations done on 8/6/19

Protocol:

Colony PCR Protocol

20 μ L Reaction

1. Prepared a PCR concentration cocktail with the following proportions: 7 μ L of dH_2O , 10 μ L PCR Mastermix, 1 μ L of the forward primer, and 1 μ L of the reverse primer.
2. Added 19 μ L of the concentration cocktail into a PCR tube.
3. Using a 10 μ L micropipette, touched the tip onto the selected colony and swirled around in the PCR tube.
4. Placed PCR tube in the thermocycler at the following generic settings:
 1. 95° C for 3:00 minutes
 2. 95° C for 1:00 minute
 3. 52° C for 1:00 minute *Annealing temperature varies depending on primer
 4. 72° C for 1:00 minute
 5. 30X (Go to Step 2)
 6. 72° C for 5:00 minutesLid Temperature: 105° C

Conclusion:

We will run the gel tomorrow.