


## Plasmid Assembly Protocol with Golden Gate Assembly

- Select a 20 nt protospacer of interest. The 3' protospacer adjacent sequence (PAM) must be NGG, where N is any nucleotide. Preference is given to:
  - Sequences with purines occupying the last four (3') bases of the protospacer.
  - Sequences on the non-coding strand.
  - Sequences in which the last 12 nt of protospacer + 3 nt PAM (15 nt total) are unique in the genome (check by BLAST with all four possible NGG sequences).
- Design two 24 nt oligonucleotides (4 nt 5' sticky end + 20 nt spacer sequence) with the sticky ends ACGC on the forward primer and AAAC on the reverse primer.
- For single spacers, anneal spacer oligos as follows:
  - Resuspend both oligos to 100 $\mu$ M in water
  - Mix 5  $\mu$ L forward primer + 5  $\mu$ L reverse primer + 90  $\mu$ L 30 mM HEPES, pH 7.8
  - Heat to 95  $^{\circ}$ C for 5 min, then ramp to 4  $^{\circ}$ C at 0.1  $^{\circ}$ C/sec
- Insert annealed spacer (or dual+spacer synthetic construct) by Golden Gate assembly.
  - X  $\mu$ L Backbone (100 ng)
  - 0.3  $\mu$ L Insert
  - 2  $\mu$ L T4 Ligase Buffer (NEB)
  - 1  $\mu$ L T4 ligase (NEB) ○ 1  $\mu$ L BbsI (NEB)
  - Fill up to 20  $\mu$ L with H<sub>2</sub>O
- Golden Gate Program:

	$^{\circ}$ C	min
9x	37.0	10
	16.0	10
	50.0	5
	65.0	20
	4.0	PAUSE



- Transform 3  $\mu\text{L}$  of each reaction to *E. coli* DH5 $\alpha$  by heat shock.
- Plate 10% of recovery culture on selective plates with 10  $\mu\text{L}$  of 0.5 M IPTG and 40  $\mu\text{L}$  of 20 mg mL<sup>-1</sup> Blue-gal (in DMSO).
- Pick white colonies to selective LB and recover plasmid.

From: iGEM Bielefeld-CeBiTec