

## Protocol for Plasmid DNA Extraction using Alkaline Lysis Method

Protocol code: Alk\_Lys

### Materials

- RNAase
- Isopropanol
- Absolute Ethanol
- Glucose
- EDTA
- NaOH
- SDS
- Potassium acetate
- Glacial acetic acid
- Tris
- Concentrated HCl
- 2mL Centrifuge tubes

### Solutions

1. **Resuspension Solution** [50 mM glucose, 10 mM EDTA(pH 8.0), 25 mM Tris (pH 8.0)]

a. Prepare the solution as described in the following table.

Reagent	Volume (mL)
1M Glucose (filter sterilized)	5
1M Tris (pH 8.0)	2.5
0.5M EDTA (pH 8.0)	1
De-ion water	91.5

- b. Autoclave or filter sterilize using membrane milipore (0.20µm).
- c. Store at 4°C.

2. **Lysis Solution** (0.2 M NaOH 1% SDS)

a. Prepare the solution as described in the following table.

Reagent	Volume (mL)
0.2M NaOH	0.2
1% SDS	1
De-ion water	8.8

3. **Neutralizing Solution** [3 M Potassium acetate (pH 6.0), 1.2 M acetic acid]

a. Prepare the solution as described in the following table.

Reagent	Volume (mL)
5M Potassium acetate (pH 6.0)	60
Glacial acetic acid	11.5
De-ion water	28.5

### Procedure

1. Grow 4 mL bacterial culture in LB medium with appropriate antibiotics overnight with shaking.
2. Transfer culture to a 2 mL tube and spin down cell culture at highest speed for 1 min at table-top centrifuge.
3. Discard the supernatant. To remove the liquid completely, place the tube upside down onto a piece of paper towel for a few seconds.
4. Add 250  $\mu$ L of Resuspension Solution into each tube and vortex to completely resuspend cell pellet.
5. Add 250  $\mu$ L of Lysis Solution and mix by gently inverting the tube 5-6 times. The solution should quickly turn colorless and become more viscous which indicates bacterial lysis has taken place.
6. Add 350  $\mu$ L of Neutralizing Solution and mix by inverting the tubes several times. At this point bacterial chromosomal DNA is usually seen as a white precipitate.
7. Centrifuge the tubes at highest speed for 10 minutes.
8. Carefully transfer the supernatant (try to not disturb the white precipitate) to a new labeled 2mL tube.
9. Add 1 volume of cold isopropanol (stored at -20 °C) to each tube and mix by inverting the tubes a few times. Spin down plasmid DNA precipitate (transparency pellet) at highest speed for 5 min.
10. Discard the supernatant and remove the remaining liquid as much as possible by leaving the tube upside-down on a piece of paper towel
11. Add 500  $\mu$ L of cold ethanol 70% (stored at -20 °C) to each tube and mix by inverting the tubes a few times. Spin down plasmid DNA precipitate (transparency pellet) at highest speed for 2 min.
12. Discard the supernatant and remove the remaining liquid as much as possible by leaving the tube upside-down on a piece of paper towel. To dry faster, keep tubes at 42°C heat block for 20 minutes.
13. Resuspend the DNA pellet with 50  $\mu$ L TE buffer + 1  $\mu$ L ARNase. Completely dissolve the pellet by pipetting solution several times.
14. Incubate at 37 °C for 30 mins.
15. Store at 4°C.

### References

Birnboim, H. C. and Doly, J. (1979). A rapid alkaline extraction procedure for screening recombinant plasmid DNA. *Nucleic Acids Res*, 7(6): 1513-1523.



Birnboim, H. C. (1983). A rapid alkaline extraction method for the isolation of plasmid DNA.  
*Methods Enzymol*, 100: 243-255.