

DATE: 19/8

LUCIFERASE CHARACTERISATION PROJECT

Plates:

Since the backbone have Chloramphenicol marker we will try to make Chloramphenicol plates. '

Recipe from iGEM

1. Recommended conc.:
Chloramphenicol 25 µg/mL
2. Add agar powder 16g/L -> autoclave
3. Add the appropriate amount of
desired antibiotic to autoclaved LB-
agar media and swirl to mix.
4. Pour ~20mL per 10cm polystyrene
Petri dish.
5. Place the lids on the plates and
allow them to cool for 30-60 minutes
(until solidified), then invert the
plates. Let sit for several more hours
or overnight.
6. Label the bottom of plates with
antibiotic and date and store in
plastic bags or sealed with parafilm
at 4°C.

400 ml of LB media ~ 20 plates

**25 µg/mL => we need 10 000 µg = 0.01 g
of Chloramph.**

PLAN

Make plates -> DNA solution from iGEM HQ
-> E-coli transformation tomorrow -> Plate
on selection media -> incubate overnight ->
Colonies containing the plasmid -> liquid
culture -> bio lector on Saturday 4 pm

DATE: 26/8

Arabinose - DH5 α experiments

DH5 α cells, transformed with biobricks (luciferase) were inoculated in LB (2 ml). Arabinose was added in different amounts. Later, CaCl₂ was added, hoping that cells might get porous by it.

Sample	Arabinose	CaCL ₂
1		0
2		15 mM
3		100 mM

These cells were later heat shocked as well. 80 C for ~2 min. But no luminescence activity recorded.
Conclusion - Should not proceed with this experiment.

DATE: 27/8

Inoculation of Top 10

We got our plate of Top 10 cells from sysbio. Single colony was inoculated in LB (5 ml)

Date: **28/8**

Top 10 cells inoculations

10 ml LB was contaminated with Top10 cells from overnight culture

Preprep for competent cells

Autoclaved tips and bottles for competent cells. 1 ml tips in special autoclave bags

SOB medium

SOB medium was prepared by the recipe. pH adjustment: 0.018 gm + 0.445 ml HCl (12 M). Sterile filtered. 250 μ l MgSO₄ and 250 μ l MgCL₂ added to 50 ml SOB medium. 100 μ l and 200 μ l Top 10 cells were added to 2 different flasks

DATE: 9/9

Responsible person: Dharmik

Luciferase expression experiment

Step 1: Grow a colony (Top 10 transformed) in 10 ml (LB+ chloramphenicol) 37°C- 12 hours.

The cells are grown in a shake flask that is kept in a beaker with ice/cold water to hold the temperature down.

Step 2: Transfer the cells to a new medium.

Contaminate 10 ml (LB+ chloramphenicol) with previously grown cells.

Grown in room temp 20°C til OD (660nm) =0,6.

Step 3: Add arabinose so that the final concentration is 0,25%

No glowing cell could be seen.

DATE: 10/9

Luciferase project

The cells are taken from room temperature shaker (was kept in ice bath overnight)

The OD was checked for the culture and measured to OD=0,9. 350µl of arabinose was added and the celles was kept in the ice bath, room temperature shaker for 30 min.

The cells were checked in the dark room and there were no glowing cells. Kept in shaker for another 30 min and still no glowing cells.

DATE: **15/9**

Luciferase Experiment

20 ml culture (LB+ chloramphenicol) +
20 μ l (25 mg/ml) chloramphenicol @ 20 c
@ 200 rpm

DATE: 17/9

Responsible person: Dharmik

Miniprep of luciferase plasmid.

Luciferase cells were inoculated on the 15th of september at 5pm but they didn't grow anything in 20°C until the 16th of september at 8pm.

The flask was moved to the 37°C incubator and kept there until 17th of September at 1 pm. *Plan: after miniprep, is it still trustworthy? Test 1, try transformation. Test 2, restriction digestion.*

The plasmid miniprep protocol was followed and the concentrations was measured to 232,4 ng/μl.

Inoculation of *E.coli* cells

E.coli cells was inoculated in LB+ chloramphenicol media, both Top10 cells and Top10+ Luciferase.

DATE: 18/9

Responsible person: Dharmik

Luciferase cells experiment

The luciferase cells are growing slower but they are growing. The non-transformed Top10 cells didn't grow in LB+ chloramphenicol.

500 µl of luciferase cells were added to 10 ml of LB + chloramphenicol at 12pm.

DATE: 19/9

Luciferase experiment

100 µl of cells from 18/9 were incubated with 20ml of fresh LB+ chloramphenicol in the shaker incubator 20°C, 180 rpm at 9 am.

No growth was seen after 12,5 h.

Responsible person: Emma

DATE: 20/9

Responsible person: Dharmik

Luciferase incubation

The growth of cells where measured on 8 am and the cells had grown. *Note: OD was not measured.* Added 1-1,5ml of 10% arabinose in 20 ml of culture. *Note: Not completely sterile, should be fine since the cells are only supposed to grow for a couple of hours.*

The working scheme over the last couple of days is summarized below:

17th Sept: Overnight culture of cells in 1 ml media

18th Sept: Overnight culture of cells in 10 ml media

19th Sept: Overnight culture of cells in

20ml media at 20 °C shaker, 180 rpm

20th Sept: added arabinose 1-1,5 ml arabinose

-> result


Restriction digestion of luciferase plasmid.

Plan: See if the biobrick is correct

- *EcoR1-PstI combo was chosen*
- *Verified in benchling, cut sites on each side of the biobrick.*

PSB1C3 backbone, the entire plasmid is 9,3kb.

Kept on 37°C for 30 min and then enzyme inactivation at 87°C for 5 min



DNA	50µl
water	30µl
Green buffer	10µl
Pst1	5µl
E.coR1	5µl

DATE: 24/9

**Luciferase digestion project:
Restriction Digestion verification**

As suggested by Oliver, All 3 enzymes in one digestion

Luciferase Biobrick - 4.5 μ l (234 ng/ μ l)

Pst1 - 1 μ l

EcoR1 - 1 μ l


Nde 1 - 1 μ l

FD Green Buffer - 3 μ l

Water - 19.5 μ l

Mixture kept @ 37 C for 15 min

Prestain Gel (1.5 μ l Gel Red), small gel,
140 V, 20 min



DATE: 29/9

Characterization experiments

Biobricks used

High Constitutive Expression Cassette (HECE)	1G, Plate 1
Orange Fluorescent Protein (OFP)	9I, Plate 4
Cyan Fluorescent Protein (CFP)	17 J, plate 6
Blue Chromoprotein (Blue)	10B, Plate 5
Terminator	E.Coli plates
GFP with degradation tag (GFP)	E.Coli plates

Biobricks in iGEM plate - Add 10 µl water, take in another tube, use 2 µl for transformation

Biobricks in Chloramphenicol plates - Inoculate cells in 10 ml chloramphenicol LB for miniprep

Transformation - 50 µl cells transformed

Date: 1/10

Miniprep

OFP	146.1 ng/μl
BFP	148.1 ng/μl
CFP	129.9 ng/μl
GFP	146.1 ng/μl
HECE	133.7 ng/μl

Biobrick - Double terminator - Plate 1, 3 D
Transformed in DH5 alpha E.Coli,
chloramphenicol medium

Make Chloramphenicol plates

12.5 mg/ml medium - 24 plates


PVU II digestion of luciferase plasmid

PVUII - 2μl

luciferase plasmid - 6 μl

FD green buffer - 2 μl

Water - 10 μl



Gel run - Very good band pattern (worth giving for sequencing)

Biobrick extraction of Terminator from iGEM plate

Date: 2/10

Extraction of Biobricks: Terminator 2

For some reason, double terminator transformation was not successful. So, retransformed it and had a backup plan. Extracted another terminator

Single terminator: 1D, plate 1
(Heads up: 10 μ l water used to resuspend DNA in plate. 2 μ l used for transformation)

Sequencing: Luciferase plasmid

15 μ l DNA + 30 μ l water: given for sequencing

Date: 7/10

Restriction Digestion of Biobricks

HCEC - 9.4 µl DNA + 1 µl EcoR1 + 1 µl
Spe1

OFP - 5.1 µl DNA + 1 µl Xba1 + 1 µl Spe1

GFP - 6 µl DNA + 1 µl Xba1 + 1 µl Spe1

CFP - 6.7 µl DNA + 1 µl Xba1 + 1 µl Spe1

Blue - 6 µl DNA + 1 µl Xba1 + 1 µl Spe1

Single Terminator - 27.4 µl + 4 µl F.D green
buffer + 2 µl Xba1 + 2 µl Pst1

T4 DNA ligase:

Use upto 5 µl of mixture in competent cells
in mixture

Buffer: 2 µl

Ligase: 1 unit

Water + DNA; 17 µl

Date: 10/10

Transformation of ligated products


Ligated products -> 5 µl in 50 µl. Rest in - 20
C

DATE: 11/10

GFP - 4 colonies, Blue chromoprotein- 1,
CFP - 0, OFP - 0

Gel Run of digests of Biobricks

- 140 V, 25 min, 1 μ l ladder, 5 μ l sample +
1 μ l FD green buffer, pre-stain (1.5 μ l Gel
Red)



GFP colonies were restreaked on new plate
and were inoculated overnight for
Fluorescence microscope view with Linnea

DATE: 12/10

Transformation of ligated products

Ligated products were re-transformed (Just to see if they did not work previously due to bad luck)

All of rest of 15 μ l ligated products in 50 μ l competent cells. Done only for CFP, GFP and Blue chromoprotein