



## Modelling – Cathepsin L activity test

### ● Aim

BioVision Cathepsin-L activity detection kit is a fluorescence-based detection technology, using AFC (7-amino-4-trifluoromethyl coumarin) labeled Cathepsin-L priority substrate sequence FR. Cell lysates or other samples containing Cathepsin-L can digest FR-AFC and release free AFC, which can be easily quantitatively detected by fluorometer or fluorescence microtiter reader. This kit can test the activity of our proteins.

### ● Materials

CL Buffer

DTT

CL Substrate Substrate Ac-FR-AFC

### ● Procedure

Based on the instructions in the Cathepsin L Activity Fluorometric Assay Kit, we designed the relevant experiments as follows:

1. Add 180ul CL Buffer into 27 wells(A1-A9, B1-B9, C1-C9) in a 96-well plate.
2. Add 0 (A1-C1), 1 (A2-C2), 2 (A3-C3), 3 (A4-C4), 4 (A5-C5), 5 (A6-C6), 6



(A7-C7), 7 (A8-C8), and 8 $\mu$ L (A9-C9) of the 10mM CL Substrate Ac-RR-AFC into plate.

3. Add 20ul cathepsin L protein into each well.

4. Read sample in a fluorometer equipped with a 400-nm excitation filter and 505-nm emission filter. Measure every hour until the fluorescence value remain unchanged.

