

Name: Rehmat Babar, Shakera Thomas

Date: 6/14/19

Goal: mini preps for pCB302 plasmid

Materials

QIAprep Spin Miniprep Kit Lot 160021667

Protocol

- a. Centrifuged the overnight culture of bacterial overnight culture at 8,000 rpm for 3 minutes at room temperature.
- b. Discarded the supernatant and resuspended the pelleted bacterial cells in 500 μ L Buffer P1 and transferred it to two different Eppendorf tube of 250 μ L each.
- c. Added 250 μ L of Buffer P2 and inverted 5 times.
- d. Added 350 μ L of Buffer N3 and immediately mixed by inverting 5 times.
- e. Centrifuged for 10 minutes at 13,000 rpm.
- f. Micropipetted 800 μ L of the clear supernatant into a spin column and centrifuged for 60 seconds and discarded the excess liquid.
- g. Added 500 μ L of PB and centrifuged the spin columns for 60 seconds. Discarded the flow through.
- h. Added 750 μ L of PE to the spin columns, centrifuged for 60 seconds, and discarded the flow through.
- i. Centrifuged the spin columns again for 60 seconds to remove residual wash buffer and discarded the flow through.
- j. Transferred the spin columns to a clean eppendorf tube and added 50 μ L of EB to the center of the spin column to elute the DNA.
- k. Allowed the spin column to stand for one minute and then centrifuged for one minute.
- l. Recorded the concentrations for each sample.

Results

Mini Prep A1 pCB302 plasmid	15 ng/ μ L
Mini Prep A1 pCB302 plasmid	5 ng/ μ L
Mini Prep A2 pCB302 plasmid	5 ng/ μ L
Mini Prep A2 pCB302 plasmid	7.5 ng/ μ L

Mini Prep A3 pCB302 plasmid	Too low
Mini Prep A3 pCB302 plasmid	7.5 ng/μL
Mini Prep A4 pCB302 plasmid	5 ng/μL
Mini Prep A4 pCB302 plasmid	2.5 ng/μL
Mini Prep A5 pCB302 plasmid	10 ng/μL
Mini Prep A5 pCB302 plasmid	5 ng/μL
Mini Prep B1 pCB302 plasmid	7.5 ng/μL
Mini Prep B1 pCB302 plasmid	7.5 ng/μL
Mini Prep B2 pCB302 plasmid	10 ng/μL
Mini Prep B2 pCB302 plasmid	10 ng/μL
Mini Prep B3 pCB302 plasmid	10 ng/μL
Mini Prep B3 pCB302 plasmid	10 ng/μL
Mini Prep B4 pCB302 plasmid	7.5 ng/μL
Mini Prep B4 pCB302 plasmid	12.15 ng/μL
Mini Prep B5 pCB302 plasmid	7.5 ng/μL
Mini Prep B5 pCB302 plasmid	10 ng/μL

Conclusion

The concentrations are too low, it is not likely the our desired plasmid was successfully transformed and/or isolated.

-

Name: Amirah Hurst

Date: 6/14/19

Goal: gel extract for Bluechromoprotien (K592009) and RFP promotor (J23102) from analytical digests performed on 6/13/19

Materials

Specific brand name and lot number

Protocol <https://www.mcdb.ucla.edu/Research/Banerjee/protocols/gelextraction--Qiagen.pdf>

1. Cut gel fragments from yesterday into 3 small pieces (1 per gel lane)
 - a. Only used 1 lane. The rest were placed in the fridge.
2. Weighed 1.5 ml microcentrifuge tube, zeroed it, placed fragment in microcentrifuge tube to get weight of fragment in mg.
 - a. K: 335.6 mg
 - b. J: 291.2 mg
3. Multiplied that weight x3 to determine ul of Buffer QG to add to the 1.5 ml microcentrifuge tube
4. Added Buffer QG
 - a. K: 1,006.8 ul
 - b. J: 873.6 ul
5. Added 30 ul QIEX 2
6. Placed tubes in 50 degrees C water bath for 10 mins vortexing every 2-3 mins to completely dissolve
 - a. Incubated for additional 5 mins because it was not completely dissolved
7. Added 1 gel volume of 2-propanol to sample, inverted to mix
 - a. K: 335 ul
 - b. J: 291.2 ul
8. Place samples in quick spin column
 - a. Used 2 spin columns for each sample
9. Centrifuge for 1 min
 - a. To bind DNA to column
10. Discarded flow through
11. Added 0.75 ml of Buffer PE to column and centrifuge for 1 min
12. Discard flow through and centrifuge for 1 min
13. Place column into a 1.5 ml microcentrifuge tube
14. Eluted DNA with 30 ul of Buffer EB, let sit for 2 mins then centrifuged for 1 min
15. Checked concentration on nanophotometer

Results:

Both samples had very low concentration

- i. J looked like maybe some was there
- ii. K was redone

Jessica and Christina did another gel extraction on blue chromoprotein (K592009) and promoter (J23102)