

Protein Extraction Protocol

Introduction

This protocol is an adaptation of a GUS (β -glucuronidase) assay provided by Zofia Dorota Jarczynska and adapted by Jacob Mejlsted

Materials

› Extraction buffer (1L)

- › 50 mM Na_3PO_4 , pH 7 (50 mL of 1 M stock solution)
- › 10 mM β -mercaptoethanol (0.7 mL of 14.4 M stock solution)
- › 10 mM Na_2EDTA (20 mL of 0.5 M Na_2EDTA stock solution)
- › 0.1% Sodium Lauryl Sarcosine (10 mL of 10% Sarcosyl)
- › 0.1% Triton X-100 (10 mL of 10% Triton)
- › MilliQ water up to 1L

› Equipment

- › Centrifuge with cooling
- › Plate reader
- › Homogenizer
- › Beads (big, metal)
- › Sterile funnels with mira cloth

› Consumables

- › Liquid nitrogen
- › Microtiterplate
- › FastPrep tubes

› Chemicals

- › Bradford reagent (like SigmaAldrich B6916)
- › BSA standard (Bovine Serum Albumin 2 mg/mL)

Procedure

Protein Extraction

1. CRITICAL Set the centrifuge at 4°C
2. Separate biomass from supernatant using mira cloth or other methods
3. Take ~100mg biomass and place in a FastPrep tube with one big metal bead
4. Place the tubes in liquid nitrogen

5. Homogenize for 1 min at 45 Hz
6. Add 500 μ l extraction buffer
7. Homogenize for 2 min at 45 Hz
8. Centrifuge for 2 min at 10,000 x g at 4°C
9. Collect the liquid phase and place in a new tube
10. Centrifuge for 15 min at 10,000 x g at 4°C
11. Aliquote the supernatant into new tubes
12. Store at -80°C

Fluorescence Assay (optional)

13. Transfer 100 μ l protein sample into a microtiter plate
14. Measure fluorescence at the appropriate wavelengths
Like (Ex/Em) 584/697 for mRFP1

Bradford Assay

15. Mix 20 μ l protein sample with 230 μ l Bradford reagent
16. Transfer to a new microtiter plate
17. Measure the blue color at 595 nm
Remember to include a standard curve and Bradford reagent
18. Determine protein concentration