

Golden gate assembly with SapI (Type IIs assembly)

Introduction

This protocol is made to assembly parts with Typell's enzyme SapI

Materials

› Reagents

- › T4 DNA ligase
- › 10X T4 DNA ligase buffer
- › SapI enzyme (NEB)
- › Reciever plasmid
- › DNA fragments with SapI overhangs
- › MilliQ water

› Materials

- › PCR tubes
- › 10uL pipette tips
- › Thermocycler

Procedure

Setting up the reaction

1. In a PCR tube, mix the following:
 - 0.5 μ L of T4 DNA Ligase
 - 2 μ L of 10X T4 DNA Ligase Buffer
 - 0.5 μ L of Type IIS restriction enzyme
 - 100 ng of receiver plasmid
 - Equimolar amounts of inserts
 - MilliQ for a total volume of 20 μ L.
2. Mix gently
3. Place the tube on a thermocycler

Thermocycler program

Table1

	A	B	C
1	Step	Temp	Time
2	1: Activation of SapI	37 °C	5 min
3	2: Activation of T4 ligase	16 °C	5 min
4	Repeat step 1 & 2 for 25 cycles		
5	3: Inactivation SapI	65 °C	20 min
6	4: Inactivation T4 ligase	85 °C	10 min
7	5: Hold	4 °C	Hold