

Name: Asma, Sijia, Jiazi, Krithika, Chiara , Laura

Date: 6/25/19

Goal:

1. Redo gel electrophoresis for PCR products
 - a. Ligations (K592009 & J23102); Vf, Vr
 - b. Pcb302 in E. Coli from papers A & B
2. Redo PCR on ligations (40 μ L reaction)
 - a. K592009 & J23102
3. Perform restriction digest on pcb302 from minipreps in E. Coli from papers A & B (colony 7) (6/20/19)
 - a. Digested with KpnI
4. Make ASP-8A medium
5. Transformation on K592009 /J23102 ligation mix from 06/17
6. O.Marina
 - a. Fed 3 mL of D. Tertiolecta

Date: 7/30/19

Goal:

1. Redo gel electrophoresis for PCR products
 - a. Ligations (K592009 & J23102); Vf, Vr
 - b. Pcb302 in E. Coli

Protocol:

Gel electrophoresis

1. Loaded 2 μ L of DNA into each well
2. Loaded 4 μ L of 1 kb plus DNA ladder
3. Ran on 2 % double comb fast E-gel for 25 minutes

Results:


Gel Key (Ligations and pcb302)

Lane #	Sample
1	Ligation 1, 150 μ L, colony 11
2	Ligation 2, 100 μ L, colony 7
3	Ligation 2, 100 μ L, colony 12
4	Ligation 1, 100 μ L, colony 9
M	1 kb plus DNA ladder
5	Pcb302, colony 7, Primers 1 & 2
6	Pcb302, colony 8, Primers 1 & 2
7	Pcb302, colony 9, Primers 1 & 2
8	Pcb302, colony 10, Primers 1 & 2

Lane #	Sample
9	Pcb302, colony 7, Primers 3 & 4
10	Pcb302, colony 8, Primers 3 & 4
11	Pcb302, colony 9, Primers 3 & 4
12	Pcb302, colony 10, Primers 3 & 4
M	1 Kb plus DNA ladder
13	Pcb302, colony 7, Primers 1 & 4
14	Pcb302, colony 8, Primers 1 & 4
15	Pcb302, colony 9, Primers 1 & 4
16	Pcb302, colony 10, Primers 1 & 4



Expected Results:



Conclusion:

This gel did not produce any solid bands.

Date: 6/25/19

Goal:

1. Redo PCR on ligations (40 μ L reaction)
 - a. K592009 & J23102

Protocol:

40 μ L PCR of Ligated BCP Samples

1. A 5x cocktail was created with the following proportions: 70 μ l diH₂O, 100 μ l PCR MasterMix, 10 μ l Vr Primer, and 10 μ l Vf Primer
2. Each PCR tube was loaded with 38 μ l of the cocktail mix and 2 μ l of the "Ligation" DNA
3. The samples were:
 - a. Ligation 1, 150 μ l, colony 11
 - b. Ligation 1, 100 μ l, colony 7
 - c. Ligation 2, 100 μ l, colony 7
 - d. Ligation 2, 100 μ l, colony 12
4. The samples were then put in the PCR Thermocycler in the Vf/Vr cycle and were left overnight

Date: 6/25/19

Goal:

1. Perform restriction digest on pcb302 from minipreps in E. Coli (colony 7) (6/20/19)
 - a. Digested with KpnI

Protocol:

30 µL Fast Digest Restriction Digest

1. Prepared a Fast Digest concentration cocktail with the following proportions: 1 µL Restriction Enzyme KpnI, 3 µL of 10X Fast Digest Buffer, and 16 µL of diH₂O.
 - a. Let enzyme thaw on ice
2. Added 20 µL of this cocktail to a clean 1.5 Eppendorf tube and then added 10 µL of DNA
3. Incubated at 37° C for 30 minutes.

Date: 6/25/19

Goal:

1. Make ASP-8A medium

Protocol:

Make ASP-8A medium

1. Vitamin 8A Mix 2X 1L

COMPONENTS	2X STOCK (1L)
p-aminobenzoic acid	0.0172g
biotin	0.001g
B12	0.001g
Choline diH2 citrate	1g
Folic acid	0.005g
Folinic acid, Ca salt	0.0004g
Inositol	2g
Nicotinic acid	0.2g
Orotic acid	0.036g
D-Pantothenic,Ca salt	0.4g
B6	0.12g
Riboflavin(B2)	0.01g
Thiamine,HCl(B1)	0.4g
Thymine	1.6g

2.ASP-8A MEDIUM 1L

COMPONENTS	1L
NaCl	25g
KCl (1M)	10ml
MgSO ₄ 7H ₂ O(1.8M)	20ml
CaCl ₂ 2H ₂ O(0.75M)	10ml
NaNO ₃ (0.58M)	1ml
KH ₂ PO ₄ (0.073M)	1ml
NTA(0.157M)	1ml
Tris Base pH 9	10ml
Adjust pH 8.5	
NH ₄ NO ₃	1ml
PII Metal Mix	10ml
8a Vitamin Mix	0.25ml
Vitamin B12	100ul
GeO ₂	2.5mg

3. Autoclaved for almost 46 minutes.

Date: 6/25/19

Goal:

1. Transformation on K592009 /J23102 ligation mix from 06/17

Protocol:

Heat Shock

1. Thawed One Shot TOP10 chemically competent cells on ice.
2. Added 2 μL of DNA sample into competent cells
3. Incubated the cells on ice for 35 minutes.
4. After the ice incubation, placed the samples into a 42° C water bath for 30 seconds.
5. **Quickly** took them out and **immediately** added 250 μL of SOC medium
6. Placed the samples into a 37° C shaking water incubator for 1 hour at 200 rpm.
7. After shaking for 1 hour, streaked 150 μL of the solution onto an agar plate with the respective antibiotics.
 - a. Ampicillin
8. Incubated plates at 37°C for at least 24 hours.