

## Transformation in *S. cerevisiae*

Works only for one plasmid transformation

- Inoculate an overnight culture of yeast strains in 10 ml YEPD and incubate at 30°C.
- Take 1.5 ml of the overnight culture and centrifuge at 4000 rpm for 3 minutes.
- Discard the supernatant and resuspend the cell pellet in 500 µl PLATE solution
- Add 2000 ng pDNA
- Incubate cells in water bath at 42°C for 30 minutes (for temperature-sensitive strains: 15 minutes incubation)
- Centrifuge at 5500 rpm for 30 s.
- Resuspend cells in 300 µl SOS medium (prepare the SOS medium shortly before use).
- Incubate cells in water bath at 30°C for 30 minutes.
- Plate cells into SD selection plate.
- Incubate the plates for a few days at 30°C (for temperature-sensitive strains: 24°C)

### SOS medium:

- 500 µl YEPD
  - 500 µl 2 M Sorbitol
  - 6.5 µl 1 M CaCl<sub>2</sub>
- For 1 ml

### PLATE:

- 20 g PEG 1500
  - 0.788 g CHC<sub>3</sub>OOLi 2 H<sub>2</sub>O
  - 0.06057 g Tris
  - 0.005 g EDTA
- For 50 ml

### YP medium:

- 10 g/l yeast extract
- 20 g/l pepton
- 20 g/l glucose

### SD-minimal medium:

- 900 ml LB
  - 100 ml C-source
  - 5 ml of each amino acid
- For 1 L

From: iGEM Bielefeld-CeBiTec