

## Online Notebook Template

Name: Laura Das Neves, Asma Khimani, Sijia, Jiazi, Amirah, Chiara Brust, Saleh Alhassam

Date: 6/18/19

Goal:

1. Simulate gel for restriction digest on Snapgene
2. Run a gel on the Restriction Digest for pCB302 samples A and B from 6/17/19 and verify if we have the correct parts.
3. Colony PCR on ligations from previous day
4. Overnight culture on ligations as well as the pCB302 plasmid
5. Create Primers for pCB302 partial sequence on snapgene

Protocol

### Colony PCR Protocol

Prepare 24 PCR tubes.

For 20  $\mu$ L Reaction

1. Prepare a PCR concentration cocktail with the following proportions: 200 $\mu$ L of diH<sub>2</sub>O, 250  $\mu$ L PCR Mastermix(2x), 25  $\mu$ L of the forward primer, and 25  $\mu$ L of the reverse primer. (Total 500uL PCR cocktail)
2. Add 20  $\mu$ L of the concentration cocktail into each PCR tube.
3. Using a 10  $\mu$ L micropipette, touch the tip onto the selected colony and swirl around in the PCR tube.
4. Place PCR tube in the thermocycler at the following generic settings:
  1. 95° C for 3:00 minutes
  2. 95° C for 1:00 minute
  3. 52° C for 1:00 minute \*Annealing temperature varies depending on primer
  4. 72° C for 1:00 minute
  5. 30X (Go to Step 2)
  6. 72° C for 5:00 minutes

Lid Temperature: 105° C

For the samples leaving in the machine, its label are as follows:

1 1:ligation 1 100  $\mu$ L colony 1,1 2:ligation 1 100  $\mu$ L colony 2,...,to 1 6  
1' 1:ligation 1 150  $\mu$ L colony1,...to 1'6  
2 1:ligation 2 100  $\mu$ L colony 1, ... to 2 6  
2'1: ligation 2 150  $\mu$ L colony 1, ... to 2'6

### Overnight Cultures



1. Add about 6 mL of LB(chloro added already) to a 15 mL Falcon tube
2. Dip a p10 tip into your selected colony and drop into the tube
3. Incubate in the water bath at 37° C at 220 rpm for 16-18 hours

## Results

S = sample

Gel 1: A samples



Lane 1	2	3	4	5	6	7	8	9	10	11	12	13	14
Gene Ruler 1 Kb plus	A1 s 1	A1s 2	A2 s 1	A2s 2	A3s 1	A3s 2	A4 s1	A 4s 2	A5s1	A5s2	empty	empty	Ladder



(run 20 min. longer)

Gel 2: B samples

Lane 1	2	3	4	5	6	7	8	9	10	11	12	13	14
Gene Ruler 1 Kb plus	empty	B1 s1	A1s2	A2s1	A2s2	A3s1	A3s2	A4s1	A4s2	B5s1	B5s2	Ladder	empty



(Run 20 min. longer)

Expected Results

