

## Protocol for Colony PCR

Protocol code: Col\_PCR

### Materials:

- DNA
- 5X Green GoTaq® Reaction Buffer
- GoTaq® DNA Polymerase (5U/ $\mu$ L)
- dNTPs
- Primers
- Nuclease Free Water
- 0.2mL tubes
- Centrifuge
- Agar plates

### Procedure:

1. Take the necessary PCR tubes for all colonies to amplify.
2. Add 25  $\mu$ L of nuclease free water to each PCR tube.
3. Resuspend the colonies in the different PCR tubes.
4. Label a petri plate with the colonies chosen and inoculate 1.25  $\mu$ L of each colony.
5. Put the PCR tubes at 95°C for 6 minutes to get the bacterium lysed and keep the plasmid in the supernatant.
6. Centrifuge at 11000 rpm for 2 minutes and 4°C.
7. Do a PCR using 2  $\mu$ L of supernatant like the DNA to amplify.
8. Continue the reacción following the protocol codified as "PCR".