

pMCY Enzyme characterization


Project: iGEM UANL 2019 Shared Project

Authors: J. Claudio Moreno-Rocha

TUESDAY, 7/23/2019

pMCY30 and pMCY90 DNA miniPrep:


Miniprepstubes.jpeg



LB media plus antibiotic

WEDNESDAY, 7/24/2019

Gel1 9030Jul21b.jpg





GEL1

For pMCY30 (BMC proteins) can chartarizate with:

*Pvu*I (3 fragment) 4.2 1.4 0.9 Kb



*Rsa*I (8 fragment) 2.3 1.4 1.1 0.6 0.5 0.3 0.09 0.05





For pMCY90 (mCherry) can chartarizate with:

PvuII (2 fragment) **2.4 0.6** Kb
RsaI (6 fragment) **1.67 0.63** 0.25 (x2) 0.18 0.048



EzDig pMCY90 (mCherry):

- PvuII Buffer B2 (Jena BioSciences)
- RsaI Buffer B2 (Jena BioSciences)

Table1				
		uL	Rx	SubTotal
1	DNA	5.0		15.0
2	Buffer	1.0		3.0
3	Enz Rsal/PvuII	0.5	X3	1.5
4	H2O	3.5		10.5
5		10.0		30.0

37oC 1-2 hrs

Gel EzDig pMCY90 (mCherry):

PvuI (30a, 30b, 30c)

RsaI (30a, 30b, 30c) (aprox 2hrs)


/---***PvuII***---/

/---***RsaI***----/

pMCY90a SE 90a 90b 90c 90a 90b 90c Mark

Agarose 0-8%, Buffer SB1X long gel 190V

Gel2 PvuII RsaI Jul21.jpg



Gel2 (500 bp DNA ladder Jena BioSciences)

THURSDAY, 7/25/2019

pMCY30 Clones

6 Clones from transformation ---> LB media + Amp (50 ng/uL) 37oC 18 hrs 900 rpm thermomixer
No bacteria

Transformation pMCY30 (Filter Paper) plus re-Growth the pMCY30b (Gel1)


Growth **30.1** and **30.3** clones in **600 uL LB** and **0.3 uL Amp** (50 ug/uL)

5 uL H₂O mQ to filter paper ---> **50 uL** Ca++ Cells to transformation ---> Lb Amp plate 37°C 18 hrs.

FRIDAY, 7/26/2019

Transformation pMCY30 (Filter Paper)
No bacteria

Re-Growth the pMCY30b




(Clone 1, 2, 3, 4)

Growth **two colonies from each clone** in 600 uL LB and without Amp

SATURDAY, 7/27/2019

Miniprep from clones pMCY30 (BMC) Clone 1, 2, 3, 4 (from bacteria culture without Amp)

pMCY90b Unk1 Unk2 30.1 30.3 31 32 33A 33B 34A 34B 3.01 3.03
Agarose 0.8% SB 1X short gel



Gel 3

Clone pMCY90 characterization *Pvu*l + *Rsa*l

*Pvu*l Buffer U (Jena BioSciences)


*Rsa*l Buffer B2 (Jena BioSciences)

Table2				
		uL	Rx	SubTotal
1	DNA	3.0		6.0
2	Buffer	1.0		2.0
3	Enz <i>Rsa</i> l/ <i>Pvu</i> l	0.5	X2	1.0
4	H2O	5.5		11.0
5		10.0		20.0

□□□/--*Pvu*ll--/ /--*Rsa*l--/

30.1 SE 30.1 32 30.1 32 Mark


Agarose 0.8% SB 1X short gel



Gel4 (500 bp DNA ladder Jena BioSciences)

MONDAY, 7/29/2019

pduU modified Synthetic DNA and promoter p3B5C in synthesis (IDT)

 Screenshot 2019-07-30 at 13.38.35.png

Estimado/a UANL iGEM 2019,

Este es un mensaje automático de Integrated DNA Technologies para informarte que tu pedido ha cambiado de estado en el Uniparts-Mexico Portal de Integrated DNA Technologies. Los detalles de tu pedido se encuentran abajo:

Fecha:7/29/2019 5:43:26 PM

Estado Actual:Approved

Observaciones del Pedido:

Nombre de la secuencia: pduUmod

Productos:gBlocks® Gene Fragments 501-750 bp

Purificación:

Secuencia:GGA TCG AAT TCG GAT TTC TAG AGG GCA TTG TAG ATA CGC TTT CGT GTT AAG AGG CTA GCG AGA ATA TTT GCA CTC CAT CCG AGA TAA GGA ATT TAC CCA TGA AAC AGG TCA CAT TGG CGC ATC TGA TTG CCA ATC CGG GTA AGG ATT TAT TTA AGA AGT TAG GCT TGC AGG ACG CCG TTA GCG CCA TCG GTA TCC TGA CTA TTA CTC CTT CTG AAG CCT CAA TCA TTG CCT GTG ACA TTG CAA CTA AAT CTG GGG CCG TAG AGA TCG GGT TCT TGG ATC GCT TTA CAG GTG CTG TAG TAT TGA CAG GTG ATG TCT CGG CAG TTG AGT ATG CCC TTA AAC AGG TCA CTC GCA CCT TAG GCG AAA TGA TGC AAT TCA CCA CAT GCT CGA TCA CAC GCA CAT AAG AGG CTA GCG AGG GTT GAC CCC AAG GGC GAC ACC CCC TAA TTA GCC CGG GCG AAA GGC CCA GTC TTT CGA CTG AGC CTT TCG TTT TAT TTG ATG CCT GGC AGT TCC CTA CTC TCG GGA CTA GTG GCT GCA GGG

Cantidad:1

Nombre de la secuencia: p3B5C

Productos:gBlocks® Gene Fragments 501-750 bp

Purificación:


Secuencia:GGA TCG AAT TCG GAT TTC TAG AGG TGT CAC CTT CAG CAA CTT TTG TTC GAT TAT CGA ACA AAT TAT TGA AAT ATC GAA CAA AAC CTC TAA ACT ACT GTG GCA CTG AAT CAA AAA ATT ATA AAC CCT GAT CAG ACA TTA CTC ATA AAA ATA CGA GGA GAG CAA CGA TGC AAC AAG AAG CAC TAG GAA TGG TAG AAA CCC GTC CTC GAA GTT CAT CAC GCG CTC CCA CTT GAA GCC CTC GGG GAA GGA CAG CTT CAA GTA GTC GGG GAT GTC GGC GGG GTG CTT CAC GTA GGC CTT GGA GCC GTA CAT GAA CTG ATT CTA CAT GAT GTA TTC GTC TAA GTG CAA TGT TAC AAT CTT AGT CAC CTT CAG CTT GTA TGT CTG GGT GCC CTC GTA GTT ACG GCC CTC GCC CTC GCC CTC GAT CTC GAA CTC GTG GCC GTT CAC GGA GCC CTC CAT GTG CAC CTT GAA GCG CAT GAA CTC CTT GAT GAT GGC CAT GTT ATC CTG GTG TTA TGA GCC ATA TTC AGG GGA CTA GTG GCT GCA GTG ATT TAG


Cantidad:1

TUESDAY, 7/30/2019

Synthetic DNA Enzimatic Pathway

MetPathway.png





Synthetic Parts:

Table3				
	A	B	C	D
1	Name	Size pb	Prefixx	Subfixx
2	Pi1VanRPcaU	1,782	EcoRI/XbaI	SpeI/PstI
3	Pi2LigM	1,617	EcoRI/XbaI	SpeI/PstI
4	Pi3AroY	1,759	EcoRI/XbaI	SpeI/PstI
5	Pi4catA	1,151	EcoRI/XbaI	SpeI/PstI
6	Pi5LigMpep	1,593	EcoRI/XbaI	SpeI/PstI

Ca++ Competent Cells *DH5α*

Overnight culture (aprox 1 mL) --> 500 mL flask with 100 mL LB media (37°C 150 rpm)

50 mL Bacteria culture --> 50 mL conic tube (2x) --> 8,000 rpm 8 mins

Supernatant out, Bacteria pellet + CaCl₂ 0.1 M --> Ice 25 mins --> 8,000 rpm 8 mins

--> Pellet + 2 mL CaCl 0.1 M ready to use

Transformation:

- pSB1C3 0.1 µg (100 ng)
- 2 µL cjBlue
- 2 µL eforRed



DNA digestion plasmid pMCY30 (BMCs) (Jena BioSciences Enzyme and Buffer)


*Xba*I Buffer Universal

*Spe*I Buffer Universal

*Pst*I Buffer Universal


*Xba*I/*Spe*I Buffer Universal

Screen Shot 2019-08-06 at 3.39.57 PM.png



XbaI SpeI X/S PstI

30.1 SE 30.1 30.1 30.1 30.1 Mark
Agarose 1.0% TBE 1X



Gel4 (500 bp DNA ladder Jena BioSciences)


WEDNESDAY, 8/7/2019

Transformation **pSB1C3 RFP** in *E. coli* DH5α Ca++ Competents Cells
37°C --> 18 hrs




Transformation **pSB1C3 *cjBlue***, **pSB1C3 *eforRed*** and **pSB3K3 RFP** in *E. coli* DH5α Ca⁺⁺
Competents Cells 37°C --> 18 hrs.

PlatesAug8.jpeg



THURSDAY, 8/8/2019


Tubes **pSB1C3 RFP** 6 clones



6 Gel clones miniprep

30.1 SE R1 R2 R4 R5 R6

Agarose 0.8% SB 1X



FRIDAY, 8/9/2019

pSB3K3 RFP miniprep 6 clones pSB1C3, pSB1C3 cjBlue 4 clones

Clones R1 R2 R4 R5 of pSB1C3 RFP characterization with EcoRI + PstI

EcoRI/PstI Buffer U (Jena BioSciences)


Table4				
		uL	Rx	SubTotal
1	DNA	1.0		4.0
2	Buffer U	1.0		4.0
3	Enz EcoRI	0.3	X4	1.2
4	Enz PstI	0.3		1.2
5	H2O	7.4		29.6
6		10.0		40.0

EcoRI/PstI

/-- pSB1C3 cjBlue --/ /-- pSB1C3 3K3 --/ /-- pSB 1C3 -/

Blue1 Blue2 Blue3 Blue4 SB1-1 SB1-2 SB6-1 SB6-2 SB6-3 SB6-4 R1 R2 R4 R5 MrK

Agarose 0.8% SB 1X 37°C 2 hrs



Gel7

MONDAY, 8/12/2019


Re-growth pMCY30 (BMC) plasmids, 6 tubes (600 µL LB + 0.6 µL Cloramphenicol, 34 mg/mL) --> 37°C-18 hrs

TUESDAY, 8/13/2019

Miniprep for 6 clones same plasmid pMCY30 (BMC)

pSB3K3 pSB1C3

Agarose 0.8% SB 1X




Gel8

Clones pSB1C3 RFP Enzyme digestion with *EcoRI* + *PstI*

EcoRI/PstI Buffer U (Jena BioSciences) 37°C 2 hrs

Table5


		uL
1	DNA	4.0
2	Buffer U	4.0
3	Enz <i>EcoRI</i>	1.0
4	Enz <i>PstI</i>	1.0
5	H2O	30.0
6		40.0



EcoRI/PstI

pSB1C3 SE (2μL) pSB1C3 (5 μL) MkR (3 μL)

Agarose 0.8% SB 1X



WEDNESDAY, 8/14/2019

pSB1C3 RFP *EcoRI/PstI* Digestion concentration in NanoDrop

485 ng/μL

Table6

	A	B	C
1		Direct Sam	1:10 Dil
2	Abs	5.815	0.683
3	260/280	1.88	1.88
4	ng/uL	391.4	48.5

THURSDAY, 8/15/2019

Pippet tips, LB Media and Amp / Cloramphenicol LB plates (30 / 34)

Get Rhamnose (0.5 gr) for encapsulation-Pep-mCherry clone of the pMCY90 plasmid.

Screenshot 2019-08-15 at 16.54.46.png

Different BMC systems alter protein expression levels and toxicity of E

To compare E expression, the different *E. coli* strains were grown in LB medium and induced under one of the following conditions: IPTG only (BMC), rhamnose only (E), or rhamnose and IPTG co-induction (E + BMC). IPTG was added to 0.5 mM at the start of growth, while rhamnose was added to 0.1 mM once cells reached an OD₆₀₀ of 0.4, about 3 h after the start of growth. The ..

Yung et al. Microb Cell Fact (2017) 16:71-88

Transformation cjBlue (13D plate 6, 2019 Kit) eforRed (11P plate 6, 2019 Kit)
DH5α Ca++ Competent Cell -- LB Agar Plates (New)

eforRedcjBlue2.jpg

No clones

MONDAY, 8/19/2019

Enzyme pathway genes and regulator

Synthetic DNA Disolution, Restriction Enzyme and ligation with pSB1C3 ligation

10 μL of H_2O to each Synthetic DNA (1 μgr per Synthetic fragment) --> **100 ngr/ μL**

Digestion with *EcoRI/PstI*

Table7						
	A	B	C	D	E	F
1	DNA	2.5		12.5		
2	Buff U 10X	1.0		5.0		
3	EcoRI	0.5	X5	2.5		
4	PstI	0.5		2.5		
5	H2O	5.5		27.5		
6		10.0		50.0	Final Con.	<u>25 ng/μL</u>

37°C/2 hrs

Ca++ Competent cell preparation

Starter culture of LB (no antibiotics) --> Grow at 37°C in shaker, overnight

(20 μL of DH5 α cells in 800 μL LB media X3)

Synthetic DNA concentration **25 ng/ μL**

Vector 485 ng/ μL --> 1:10 dilution --> **48.5 ng/ μL** (2,033 pb)

Table8						
	A	B	C	D	E	F
1		Pi1VanRPcaU	Pi2LigM	Pi3AroY	Pi4catA	Pi5LigMpep
2	2,033 pb Vec	1,782 pb	1,617 pb	1,759 pb	1,151 pb	1,593 pb
3	DNA vector	1.00	1.00	1.00	1.00	1.00
4	DNA Insert	5.11	4.64	5.04	3.30	4.57
5	Buffer Lig 10X	2.00	2.00	2.00	2.00	2.00
6	Enzyme T4 Lig	0.25	0.25	0.25	0.25	0.25
7	H2O	11.64	12.11	11.71	13.45	12.18
8		20.00	20.00	20.00	20.00	20.00
9	<i>Vector</i>	<i>48.5 ng</i>	<i>48.5 ng</i>	<i>48.5 ng</i>	<i>48.5 ng</i>	<i>48.5 ng</i>
10	<i>Insert</i>	<i>127.8 ng</i>	<i>116.0 ng</i>	<i>126.1 ng</i>	<i>82.5 ng</i>	<i>114.2 ng</i>

DNA ligation transformation

Table9		
	A	B
1	pSB1C3 (control)	2 uL (aprox 0.1 ng)
2	pi1Van	5 uL (44.0 ng)
3	Pi2LigM	5 uL (41.1 ng)
4	Pi3AroY	5 uL (43,6 ng)
5	Pi4catA	5 uL (32.7 ng)
6	Pi5LigM	5 uL (40.6 ng)

DH5α pre Culture (800 µL LB media no antibiotic) 37°C 900 rpm X2

WEDNESDAY, 8/21/2019

No colonies.

New Ca++ Competen Cells --> Transform pSB1C3RFP as control (aprox 5 ng)

Ca++ Competent Cells from Cellular Biology Lab --> LB Cm Plates


Table10		
	A	B
1	pSB1C3 (control)	1 uL (aprox 2 ng)
2	pi1Van	5 uL (44.0 ng)
3	Pi2LigM	5 uL (41.1 ng)
4	Pi3AroY	5 uL (43,6 ng)
5	Pi4catA	5 uL (32.7 ng)
6	Pi5LigM	5 uL (40.6 ng)

THURSDAY, 8/22/2019

Clone in plate Pt2, grow it in LB (Cm) 600 µL (0.6 µL Cm 1000X) 37°C 900 rpm (Thermomixer Eppendorf)

FRIDAY, 8/23/2019

Clone in plate Pt4, grow it in LB (Cm) 600 µL (0.6 µL Cm 1000X) 37°C 900 rpm (Thermomixer Eppendorf)



clone Pi4catA


DNA miniprep and *EcoRI/PstI* digestion (microwave oven)

Table11

	A	B	C	D
1	DNA	2.0		4.0
2	Buff U 10X	1.0		2.0
3	EcoRI	0.3	X2	0.6
4	PstI	0.3		0.6
5	H2O	6.4		12.8
6		10.0		20.0

5 seconds in a regular microwave oven

-- pSB1C3 RFP -- /EcoRI/PstI/
Contr1 Contr2 Contr3 Pi4 Ctrl Pi4 MrK
Agarose 0.8% SB 1X



Gel 10 Pi4 miniprep

SUNDAY, 8/25/2019

Pick colonies from **Pi2** (**Pi2LigM**) and **Pi4** (**Pi4catA**)
LB (Cm) 600 µL (0.6 µL Cm 1000X) **37°C 900 rpm**

MONDAY, 8/26/2019


Miniprep **Pi2** and **Pi4** clones

DNA miniprep and **EcoRI/PstI** digestion (microwave oven)

Table12				
	A	B	C	D
1	DNA	1.0		3.0
2	Buff U 10X	1.0		3.0
3	EcoRI	0.3	X3	0.9
4	PstI	0.3		0.9
5	H2O	7.4		22.2
6		10.0		30.0

10 seconds in a regular microwave oven

pSB1C3RFP /- *EcoRI/PstI* -/
Contr3 **Pi2** **Pi4** **Ctrl** **Pi2** **Pi4** MrK
Agarose 0.8% SB 1X



Pick colonies from **Pi1** (**Pi1Van**) and **Pi3** (**Pi3AroY**)
LB (Cm) 600 µL (0.6 µL Cm 1000X) **37°C** 900 rpm

TUESDAY, 8/27/2019

Miniprep **Pi1** and **Pi3** clones

pSB1C3RFP

Contr3 **Pi1** **Pi3**

Agarose 0.8% SB 1X


 Gel 12 Pi1Pi3SE 1b.png

pSB1C3 RFP

Contr3

Pi1

Pi3



Agarose 0.8% SB 0.8X

Gel 12 Clones Pi1 and Pi3

iGEM Kit [Backbones vector](#)

pSB1C3 backbone --> 25 ng/μL

125 μg Digested with *EcoRI* + *PstI* --> 37°C 2 hrs

DNA Ligation 2

Table13						
	A	B	C	D	E	F
1		Pi1VanRPcaU	Pi2LigM	Pi3AroY	Pi4catA	Pi5LigMpep
2	2,033 pb Vec	1,782 pb	1,617 pb	1,759 pb	1,151 pb	1,593 pb
3	DNA vector	2.0	2.0	2.0	2.0	2.0
4	DNA Insert	3.3	4.7	5.2	3.4	4.7
5	H2O to 5uL	1.7	0.3	0.0	1.6	0.3
6	Buffer Lig 10X	2.00	2.00	2.00	2.00	2.00
7	Enzyme T4 Lig	0.5	0.5	0.5	0.5	0.5
8	H2O	10.5	10.5	10.5	10.5	10.5
9		20.00	20.00	20.00	20.00	20.00
10	Vector	131.5 ng	119.3 ng	129.8 ng	84.9 ng	117.5 ng
11	Insert	50 ng	50 ng	50 ng	50 ng	50 ng

25°C overnight

WEDNESDAY, 8/28/2019

DNA ligation 2 transformation in DH5α

Pi1Van -- Pi2LigM -- Pi3aroY -- Pi4catA -- Pi5LigMpep -- Control (100 pg/μL)

--> LB Agar Cm 37°C overnight

Set Seed Culture DH5α, 2 Tubes 800 μL LB (No antibiotic) --> 37°C overnigh

THURSDAY, 8/29/2019

DH5α Ca++ Competent Cells transform a plasmid pSB1C3mCherry (4 ng/μL) as control

FRIDAY, 8/30/2019

No colonies, Try new Competent cells

MONDAY, 9/2/2019

Re growth clones pMCY30 (BMC) and pMCY90 (mCherry-Internalization peptide) bacteria culture in LB.

TUESDAY, 9/3/2019

pMCY30 (BMC) and pMCY90 (mCherry-Pep) growth LB media

--> LB Agar Cm 37°C overnight

Miniprep to 30.1 30.2 90₁ 90₄ RFP₁ RFP₂ (pSB1C3RFP)

Pick new clones from LB Plates:


Set 1 30₁ 30₂ 30₃ 30₄ Set 2 30₁₁ 30₁₂ 30.3 30.4

Miniprep to Set1 --> PvuII Digestion

Table14 ^

	A	B	C	D
1	DNA	1.0		4.0
2	Buff U 10X	1.0		4.0
3	PvuII	0.3	X4	1.2
4	H2O	7.7		30.8
5		10.0		40.0

/-- PvuII --/ E/P SE
30₁ 30₂ 90₄ RFP RFP




No positive digestion

Ca++ Competent Cells Test

DNAs used for test:

- a) pMCY30 (Amp)
- b) pMCY90 (Kan) DNA --> 1 : 1000 dil aprox. 4 ng/ μ L
- c) pSB1C3 RFP (Cm)

1 μ L (1 : 1000 dil) each plasmid --> 50 μ L Bacteria --> ice, 20 min --> 42°C 90 sec --> 200 μ L LB media --> 37°C shaking 20 min --> LB Plate --> 37°C overnigh.




TUESDAY, 10/8/2019

Clone pMCY90 (mCherry-Internalization Pep) induction


Overnigh growth --> 1:10 dil (2,700 μ L LB media + 300 μ L Bac)
--> 15 mL conic tube/900 rpm/room temp

Clone90Growth.jpeg




Clone pMCY90 pre-induction

pMCY90RhaTubes.jpeg




Pre-induction

pMCY90RhaTubesb.png



Induction Rha (0.5 mM) + Control

2 hrs growth --> + 3 μ L Rhamnose 500 nM (Final conc 0.5 mM) + 3 hrs 900 rpm
 --> check fluorescence (587 nm ext / 610 nm emiss)



THURSDAY, 10/10/2019

Synthetic DNA design from each circuit construction

All parts of the genetic circuit were synthesized in IDT recibed lyophilized



The individual parts were dissolved in H₂O molecular grade. They were used to make the digestions and subsequently clones in to pSB1C3 RFP vector.

SUNDAY, 10/13/2019

Clones from Pi's fragments ligation in pSB1C3 RFP vector


Clones from the constructions:

Pi2 LigM no integration peptide

Pi3 aroY + integration peptide

Pi4 catA + integration peptide

EcoRI + *PstI* digestion



Gel 14 Pi's parts digest

WEDNESDAY, 10/16/2019

Pi4 *catA* clone characterization.

Clone 4L: enzyme characterization with *PvuII* enzyme



Gel 15 Pi4 catA characterization

WEDNESDAY, 10/23/2019


catA enzyme expression in *E. coli*

In order to express the **catA** enzyme in bacteria, the construction was cloned with **three constitutive promoters** from the **Anderson** collection.


J23100 2547 au

J23102 2179 au

J23104 1831 au




Cloning strategy



Vectors digestion: J23100, J23102 and J23104 with *SpeI* and *PstI*

Table15				
	A	B	C	D
1	DNA	2.0		6.0
2	Buff U 10X	2.0		6.0
3	SpeI	0.3	X3	0.9
4	PstI	0.3		0.9
5	H2O	15.4		46.2
6		20.0		60.0




Gel 16 Promoter *SpeI/PstI*

Vectors digestion: Pi4 catA with *SpeI* and *PstI*

Table16				
	A	B	C	D
1	DNA	2.0		4.0
2	Buff U 10X	2.0		4.0
3	<i>SpeI</i>	0.3 x2		0.6
4	<i>PstI</i>	0.3		0.6
5	H2O	15.4		30.8
6		20.0		40.0





Gel 17 Pi4 catA XbaI/PstI

Pi4 catA and Constitutive promoters ligation

Table17				
	A	B	C	D
1	DNA catA	2.0		6.0
2	DNA vector	1.0		3.0
3	Buffer 10X lig	2.0	x3	6.0
4	T4 Ligase	2.0		6.0
5	H2O	13.0		39.0
6		20.0		60.0

DNA ligation transformation in Top10

J23100 + catA -- J23102 + catA -- J23104 + catA (100 ng/μL)

--> LB Agar Amp 37°C overnight

Single-temperature Double Digest

Introduction

This is the Double Digest Protocol with Standard Restriction Enzymes, using a common reaction and same incubation temperature for both enzymes.

More information from NEB can be found [here](#).

Double Digests can be designed using [NEB's Double Digest Finder](#).

See the [NEBuffer Activity/Performance Chart with Restriction Enzymes](#) for the incubation temperatures.

[NEBcloner](#) will help guide your reaction buffer selection when setting up double digests.

Materials

- › DNA 1 µg
- › NEBuffer
 - › 1X
- › NEB Restriction Enzymes
- › Deionized Water

Procedure

Single Temperature DD Reaction

- ✓
1. Set up the following reaction (total reaction volume 50 µl).

Table2			^
	A	B	
1		Reagent Volumes (µl)	
2	Buffer (10x)	5	
3	DNA *	Input Volume for ng	
4	Restriction Enzyme #1 **	1	
5	Restriction Enzyme #2 **	1	
6	Deionized Water (µl)	48	
7	Total Volume (µl)	50	

* Recommended maximum of 1 µg of substrate per 10 units of enzyme.
** Restriction Enzymes should be added to the mixture last.

- ✓ 2. Mix components by pipetting the reaction mixture up and down, or by "flicking" the reaction tube.
- ✓ 3. Quick ("touch") spin-down in a microcentrifuge. Do not vortex the reaction.
- ✓ 4. Incubate for 1 hour at the enzyme-specific appropriate temperature.

Can be decreased to 5-15 minutes by using a [Time-Saver™ Qualified Restriction Enzyme](#)

See the [NEBuffer Activity/Performance Chart with Restriction Enzymes](#) for the incubation temperatures.