

QIAprep® Spin Miniprep Kit (Plasmid Isolation)


Catalog Nos. 27104 and 27106

Centrifugation Protocol

All centrifugation steps are carried out at 13,000 rpm ($\sim 17,900 \times g$) in a conventional table-top microcentrifuge.

The following procedure should be performed at room temperature (15-30°C).

1. Pellet 1-5 ml bacterial overnight culture by centrifugation at >8000 rpm ($\sim 6800 \times g$) for 3 min at room temperature (15-25°C).
2. Resuspend pelleted bacterial cells in 250 μ l Buffer P1 and transfer to a microcentrifuge tube.
3. Add 250 μ l Buffer P2 and mix thoroughly by inverting the tube 4-6 times until the solution becomes clear. Do not allow the lysis reaction to proceed for more than 5 min. If using LyseBlue reagent, the solution will turn blue.
4. Add 350 μ l Buffer N3 and mix immediately and thoroughly by inverting the tube 4-6 times. If using LyseBlue reagent, the solution will turn colorless.
5. Centrifuge at room temperature at 13,000 rpm ($\sim 17,900 \times g$) in a table-top microcentrifuge.
6. Apply 800 μ l supernatant from the last step to the QIAprep 2.0 spin column by pipetting. Centrifuge for 30-60 s and discard the flow-through.
7. **For strains with endA+ or other bacteria strains with high nuclease activity or carbohydrate content:** Wash the QIAprep 2.0 spin column by adding 0.5 ml Buffer PB. Centrifuge for 30-60 s and discard the flow-through.
8. Wash the QIAprep 2.0 spin column by adding 0.75 ml Buffer PE. Centrifuge for 30-60 s and discard the flow-through.
9. Centrifuge for 1 min to remove residual wash buffer.



10. Place the QIAprep 2.0 column in a clean 1.5 ml microcentrifuge tube. To elute DNA, add 50 μ l Buffer EB (10 mM TrisCl, pH 8.5) or water to the center of QIAprep 2.0 spin column, let stand for 1 min, and centrifuge for 1 min.

From: iGEM Bielefeld-CeBiTec