

PCR

Introduction

PCR for the amplification of Insert-DNA.

Materials

- Polymerase
- Polymerase-Buffer
- dNTPs
- Primer
- Template
- nuclease free Water

Procedure

1. Pipette the following volumes:

	Volume [μ l]
Q5-Buffer	10
Q5-Polymerase	0.5
dNTPs, 2mM	5
Primer fw, 10μM	5
Primer rv, 10μM	5
Template-DNA (10 ng/Probe)	1
nuclease-free H2O	23.5
Complete Volume	50

2. Distribute 50 μ l of the PCR-Mix to one PCR-Tube. Amplify the DNA in the Thermocycler, use the following programme:

Step	Temperature	Time
		in s e

		c
Initial Denaturation (25-35 cycles)	98°C	1 2 0
Denaturation	98°C	2 0
Annealing	64°C	3 0
Extension	72°C	6 0
Final Extension	72°C	6 0 0
Cooling	4°C	f o r e v e r

3. Purify with pRC-CleanUp or Prep-Gel.