

Bacteriophage Enrichment

Introduction

This protocol describes the *P.larvae* bacteriophage enrichment process in order to reproduce bacteriophages with a high titer ($>10^9$ PFU/mL).

Materials

- BHI Media
 - BHI Powder 37 g/L
 - Yeast Extract 3 g/L
- 1x TM Buffer
- sterile H₂O

Procedure

- **Day 1 - *Paenibacillus larvae* (*P.larvae*) Bacterial Lawn Preparation**
 1. Take one row of the bacteria lawn from a fresh (max 2 days old) *P.larvae* streak plate and dissolve them in 200 μ L sterile dd.H₂O
 2. Plate out those 200 μ L of the bacterial suspension on a MYPGP Agar plate
 3. Incubate overnight, 37 °C
- **Day 2**
- 5 mL *P.larvae* ONC:
 1. Resuspend the *P.larvae* ONC Plate with 2 mL ddH₂O (sterile) and pipette it in a 1.5 mL reaction tube
 2. Centrifuge 5,000 x g, 5 min in the eppifuge
 3. Discard the supernatant, and dissolve the Pellet in ca. 250 μ L ddH₂O (sterile)

4. Create 1:50 and 1:500 dilutions from the Probe and measure OD_{600}
 5. Set the OD_{600} of the bacterial probe to be $OD_{600} = 15$
 6. Pipette 5 mL BHI Media and 50 μ L of bacterial suspension ($OD_{600} = 15$) into a 15 mL falcon tube
 7. Add 500 μ L Phage lysate to the falcon tube
 8. Incubate overnight at 37 °C, 200 rpm
- 50 ml *P.larvae* ONC
 1. Take one row of the bacteria lawn from a fresh (max 4 days old) *P.larvae* streak plate and dissolve them in 200 μ L sterile dd.H₂O
 2. Pipette the bacterial suspension into 50 mL BHI Media
 3. Incubate overnight at 37 °C, about 200 rpm
 - **Day 3**
 1. Add 500 μ L CHCl₃ (Chloroform) to the Probe (5 mL ONC), vortex and let it sit for 5-10 min
 2. Centrifuge the CHCl₃-free suspension (2 mL is enough) at 12,000 x g, 15 min
 3. 50 mL *P.larvae* ONC + 1 mL Phage lysate
 4. Incubate at RT, 15 min
 5. Add the Probe into 400 mL BHI Media
 6. Incubate at 37 °C, ca. 120 rpm, 4 days

- **Day 7**

1. Cool the ultracentrifuge rotor at 4 °C
 2. Give up to 0.5 M NaCl and 5-10 % CHCl₃ to the Phage culture
 3. Incubate 5-10 min, ca. 200 rpm
 4. Transfer the CHCl₃ free solution to the 25 mL ultracentrifuge tubes
 5. Centrifuge 12,000 x g for 15 min
 6. Transfer the Phage-supernatant in a sterile 500 mL Media bottle
 - Store at 4 °C
 - If not using that day, put 1 drop of CHCl₃ in it in order to prevent bacterial growth
 7. Wash the ultracentrifuge tubes with soap, and let them dry
 8. Take 1 mL Aliquot for the phagetiter determination
 9. Centrifuge the Phage-supernatant at 100,000 x g, 4 °C, 1.5 h
 10. Discard the supernatant and fill up the tubes with the rest of the Phage lysate
 11. Centrifuge and repeat at above until there is no more phage lysate
 12. Discard the supernatant
 13. Pipette about 3 mL of sterile 1x TM Buffer into the tubes
 14. Put the tubes (45° angle) into a shaker, overnight, 16 °C, 140 rpm in order to dissolve the phage pellet
- **Day 8**
 15. Dissolve the undissolved phage pellet with a inoculating loop
 - store at 4 °C
 16. Unite all the phage pellet solutions into one sterile 50 mL falcon tube

17. Determine the phage titer via Plaque drop assay
18. If the phage titer is $>10^9$ PFU/mL, than proceed to the purification step with a caesium density gradient centrifugation.