



COLOR CODING KEY

Characterization

Construct

Data Collection/Analysis

Improvement of Biobrick

Plant-Care

Plasmid

Plasmid/Construct Design

Preparatory Work

Sunday August 4th

- Resuspended agrobacterium
 - OD600 of agro 25mL culture is 1.925
 - Pelleted the agrobacterium in a centrifuge for 5000RPM for 15 minutes
 - Resuspended pellet in agrobacterium resuspension buffer in a 15mL tube but then moved to a 250mL flask
 - Later on, adjusted OD600 to 0.4
 - Added 100uL of Acetosyringone into the buffer
 - Flask placed into shaking incubator for 28C at 200RPM until agroinfiltration
- Agroinfiltration of 10 plants:
 - 3 controls with PBS and agrobacterium resuspension buffer
 - 2 infected with TMV only
 - 2 infiltrated with agrobacterium only
 - 3 infected with TMV and infiltrated with agrobacterium
 - Experimented leaves are marked with nail polish

Monday August 5th

- Gel Electrophoresis:
 - Good results: bands are ~1kb

Tuesday August 6th

- Nanodrop 7 tubes
 - Average concentration around 450ng/uL
-

Wednesday August 7th

- Liquid inoculation of the following:
 - N gene
 - Silencing suppressor
 - Hifi restreak
 - Made 7 KAN antibiotic plates
 - Resuspended gold characterization materials
 - Followed protocol from IDT
 - PCR'd gold characterization promoters
 - Used primer 67 and 72 (35s forwards and XRN1-2 reverse)
 - Followed protocol for Q5 mastermix
 - Left in thermocycler overnight
-

Thursday August 8th

- Nanodrop:
 - Nanodropped all gold characterization
 - Concentration ~400ng/uL
- RE digest (P2,T7, Backbone [-20])
 - Reagents
 - Xba1 - 1uL
 - Spe1 - 1uL
 - Cutsmart - 5uL
 - DNA - 1uL
 - Nuclease Free Water to a total of 25uL
- Ligation
 - Reagents
 - 10X T4 Ligase Buffer - 2uL
 - Backbone DNA - 1.5uL
 - Inserts (P1,P2,P3,T7) - 2uL
 - T4 Ligase - 1uL
 - Nuclease Free Water - 13.5uL
 - Same protocol
 - Used Xba1 an Spe1
- Transformation
 - Transformed 4 ligations
 - T7

- with the first two enzymes
- and the second two enzymes
- P2
 - with the first two enzymes
 - and the second two enzymes

Friday August 9th

- Lab was cleaned
- Three 1% gel electrophoresis were ran on the liquid inoculated:
 - N gene
 - Silencing suppressor
 - Hifi restreak
- Used HiFi on 3 of 4 gold characterization promoters
 - Hifi of pcr'd backbone + promoter + mCherry + tNOS
 - PCR of backbone was used directly
 - All inserts were diluted to a point where 1uL was used
 - 5x concentration was used for any piece smaller than 250 bases
- Transformation:
 - Transformed E.Coli with HiFi products and plated
 - GC 1- 35s biobrick from registry
 - GC-3 stress response promoter
 - GC-4 heat shock response promoter
 - Positive control on AMP plates

Saturday August 10th

- Agroinfiltration with the following:
 - N gene
 - Silencing suppressor
- Miniprepped agrobacterium
- Transformed agrobacterium
 - Placed in a shaking incubator at 28C for 175RPM
- Agroinfiltrated:
 - 3 Controls
 - 2 TMV
 - 2 Agrobacterium
 - 3 Agrobacterium + TMV