



COLOR CODING KEY

Characterization

Construct

Data Collection/Analysis

Improvement of Biobrick

Plant-Care

Plasmid

Plasmid/Construct Design

Preparatory Work

Sunday July 14th

- Created and autoclaved LB media
- Streaking:
 - Streaked plates with 500uL of AMP and 500mL LB broth
 - 100uL of Hifi/Positive control transformed cells
- Transformation:
 - Used HiFi kit
- Plants:
 - All inoculated leaves have similar numbers of inoculation sites visible
 - Old inoculated plants had stunted growth
 - Sheared off big leaves for visualization of stem and other leaves

Monday July 15th

- Liquid inoculation of plates from the previous day
- Plants:
 - TMV has not entered the major leaf vasculature yet

Tuesday July 16th

- Liquid cultures taken out
 - Cultures appeared clear
 - Cultures did not pellet when centrifuged
- Plates:
 - Made CAM plates

- Prepare the CAM from the powder in the -20C fridge
 - Plants:
 - Transplanted 2-week sprouts into 16 square wells and circular ports + osmocote
 - Thinned the plants and sowed another batch
 - Sub Irrigated with fertilizer
-

Wednesday July 17th

- Liquid inoculation:
 - 3 tubes of HFA
 - 3 tubes of HFB
 - 1 tube for positive control
 - 1 tube for positive control A
 - 1 tube for positive control B
 - 1 tube for negative control
 - Created YPD agar plates
 - Used 13g of YPD agar powder in 200mL MilliQ water
 - Poured 6 plates of YPD agar
 - Streaking:
 - 1 control plate
 - 2 of Rowley's XRN1 knockout plate
 - 2 XRN1 knockout plate
 - Plates were placed into the incubator at 35C
 - Plants:
 - Created our own inoculation solution
 - Gained five 2 mL tubes of freeze dried ground up leaves
-

Thursday July 18th

- Liquid inoculations:
 - None of the positive controls worked
- Miniprep the liquid inoculations
- Nanodrop:
 - Results varied for both HiFi B inoculations
- RE Digest:
 - Used 40uL nuclease free water, 4uL DNA, 5uL cutsmart, and 1uL EcoRV enzyme
 - Placed into thermocycler at 65C for 20 minutes
- Created 16 plates
- Cleaned the hood with bleach and placed under UV light for 20 minutes

- PCR
 - Used 2x Q5 mastermix protocol for 25 uL reaction
 - For template, a double digested plasmid for Prelim plasmid extension PCR was used
 - For construct plasmid, a PmeI digested plasmid was used
 - Tubes were labeled L-PVT and L-PMY
 - Resuspension:
 - Resuspended DNA
 - Created 3 promoter plates and 1 GFP plate
 - Promoters: **BBa_K895000**, **BBa_K541503**, **BBa_R0080**
 - GFP: **BBa_E0040**
 - Stored in -20C fridge
-

Friday July 19th

- Plates
 - Threw out old plates
 - Placed plates in incubator into the 4C fridge
 - Created and autoclaved KAN and CAP LB broths
- Nanodrop:
 - Nanodropped both PCR constructs (L-PVT and L-PMY)
- Gel Electrophoresis:
 - Added 400uL of TAE and 12uL of Ethidium bromide to fill gel casket
 - For loading, we used 2uL DNA, 3uL water, and 1uL of loading dye
 - Ran for 100V
 - Made two 1% gels
 - 60mL TAE, 1.8uL of Ethidium bromide and 0.6g agarose for each gel
- Opentrons Collaboration
 - Started the human protocol for the fluorescein salt
 - Did Opentrons protocol for fluorescein salt and Abs600
- Organized the 4C fridge
- RE Digest
 - RE digested the big miniprep stock for both prelim and construct
 - Prelim
 - Hind III
 - Sac 1
 - Construct
 - Pme 1
- PCR

- PCR of main project plasmid and extension PCR of pre lim plasmid
 - RE digest product used as template for PCR
 - Followed protocol for 2X Q5 master mix for 25 uL reaction
- Used program CPL-58 for main project plasmid
- Used program EXT-66 for prelim plasmid
- Left in thermocycler and held at 4C overnight
- **Streaking:**
 - 3 AMP plates
 - 2 streaked with PVT
 - 1 control streaked with KAN resistant bacteria (PMY)
 - 3 KAN plates
 - 2 streaked with PMY
 - 1 control streaked with AMP resistant bacteria (PVT)
- **Liquid inoculation**
 - 6 liquid cultures and 1 control
 - 2 EHA 105 in LB
 - 2 GV 3101 in LB + RIF (10ug/mL)
 - 2 GV 2260 in LB + RIF (10 ug/mL)
 - Placed into the incubator at 28C and 250RPM
- **Transformation:**
 - Transformed promoter 1,2,3 and GFP into DH5-alpha cells
 - Used 1uL of resuspended DNA from the plates
 - For controls
 - resuspended the DNA from the competency test kit
 - These were positive controls
- **Plates:**
 - Promoter 1,2,3 on CAM plates
 - Diluted and concentrated
 - GFP on AMP plates
 - Diluted and concentrated
 - Controls
 - Positive
 - R and O on CAM plates
 - Negative
 - R and O on AMP plates
 - GFP on CAM plates

Saturday July 20th

- **Plates:**

- PMY with AMP plates control had growth
 - PVT with KAN plates control had growth
 - Made two 70% ethanol bottles
 - Gel Electrophoresis
 - Used 400mL TAE and 12uL Ethidium bromide
 - Used a 1% gel
 - Used 3uL water, 2uL DNA, 1uL loading dye
 - Placed 4uL into each well
 - Used 100V to run gel
 - Made another three 1% gels
 - Nanodrop - good results
 - HiFi:
 - Followed HiFi protocol with 0.5 total pmols of DNA
 - PCR plasmid DNA was used directly
 - Insert DNA was diluted until 0.5uL would need to be added to the reaction
 - Ran reaction for 2 hours
 - Transformation:
 - Followed transformation protocol
 - NEB protocol recommended tenfold serial dilutions and plating those results
 - Plates labeled 0.1x were diluted before plating
 - Then placed into an incubator
 - Plates:
 - Small growth in negative controls
 - Poor growth on GFP plate
 - Liquid inoculation:
 - Promoter 1 (x and y)
 - Promoter 2 (x and y)
 - Promoter 3 (x and y)
 - GFP (x and y)
 - Controls
 - GFP in CAM media
 - Promoter 3 in AMP media
-