

Plaque Drop Assay

Introduction

With this assay you can determine the phage titer concentration

Materials

- sterile MYPGP Agar plate
- sterile dd. H₂O - 30 mL
- BHI Media - 20 mL to dissolve the phages
 - BHI Powder 37 g/L
 - Yeast Extract 3 g/L
- BHI Top agar - 20 mL
 - BHI Powder 37 g/L
 - Yeast Extract 3 g/L
 - 0.7 % Agar

Procedure

- Day 1
 1. Take one row of the bacteria lawn from a fresh (max. 2 days old) *P.larvae* streak plate and dissolve them in 200 µL sterile dd.H₂O
 2. Dissolve the bacteria (resuspend, vortex)
 3. Plate out the 200 µL bacteria suspension onto a MYPGP agar plate
 4. Incubate overnight at 37 °C
- Day 2
 1. Prepare and autoclave BHI Top Agar and cool it at 55 °C
 2. Turn the photometer on
 3. warm the MYPGP Plates at 37 °C

4. Resuspend the *P.larvae* ONC Plate with 2 mL ddH₂O (sterile) and pipette it in a 1.5 mL reaction tube
 5. Centrifuge 5,000 x g, 5 min in the eppifuge
 6. Discard the supernatant, and dissolve the Pellet in ca. 250 µL ddH₂O (sterile)
 7. Create 1:50 and 1:500 dilutions from the Probe and measure OD₆₀₀
 8. Set the OD₆₀₀ of the bacterial probe to be OD₆₀₀ = 15
 9. Pipette 50 µL OD₆₀₀ = 15 bacterial suspension into a glass epprouvette, add 4 mL BHI-Top agar (~55 °C), mix, pour it quickly on a MYPGP agar plate
 - the plate should look smooth (no air bubbles)
 10. Cool the MYPGP Plate at RT for 10 min
 11. Prepare the phage dilutions (dilute in BHI media)
 - Dilutions from 10⁻¹ to 10⁻⁹
 12. Pipette a 5 µL drop from the phage dilution on the cooled MYPGP-Top agar plate
 13. Let it sit at RT for 20-30 min
 14. Incubate overnight at 37 °C, do not invert the plate
- Day 3
 15. Count the plaques
 16. Determine the Phage titer
 - PFU/mL = (Plaque Count * 200) / Dilution