

Agrobacterium tumefaciens Mediated Transformation of pCB302-gfp-MBD into S. microadriaticum via Bead Beating

1. Placed Symbiodinium cells in 1 mL culture medium (ASP-8A) in a 2mL cryotube containing a dry volume of 200  $\mu$ L (about 500 mg) acid-washed, sterile glass beads.
  - a. Mini-bead beater [S. Micro]=  $3.0 \times 10^5$  total cells from flask: S. Micro ASP-8A 8/2/19 100 mL; 2
  - b. Lambert Bead Beater [S. Micro]=  $9.6 \times 10^5$  total cells from flask: S.M. 75 mL + 1mL ASP-8A 8/21/19
2. Then, added 350  $\mu$ L of 20% polyethylene glycol (PEG) to the suspension.
3. Shook the tube in a bead beater at 4200 rpm for 90 seconds.
  - a. Mini-bead beater [S. Micro]=  $1.6 \times 10^4$  cells/mL
  - b. Lambert Bead Beater [S. Micro]=  $6.0 \times 10^4$  cells/mL
4. After shaking, transferred the cells to a new sterile 2 mL cryotube and washed to remove the PEG by pelleting at 3,000 g for 3 minutes then resuspending in 1 mL ASP-8A.
  - a. Mini-bead beater [S. Micro]= 0 visible cells
  - b. Lambert Bead Beater [S. Micro]=  $4.4 \times 10^4$  cells/mL (some seen swimming)
5. Added 150  $\mu$ L of Agrobacterium culture (OD600= 1.5) harboring the pCB302 plasmid.
6. Incubated in fresh ASP-8A medium without antibiotics in the dark for 2 days before selection.
7. To select, added 50  $\mu$ g/mL Kanamycin and 50  $\mu$ g/mL Ampicillin. In addition, 50  $\mu$ g/mL Kanamycin was always present to prevent any bacterial growth during the selection process.
  1. 50  $\mu$ L of 1000x kanamycin
  2. 25  $\mu$ L of 50  $\mu$ g/mL Ampicillin