

Name: Chiara Brust

Date: 8/31/19

Goal:

1. Glycerol stocks of mCherry from 8/29/19 transformation

Name: Chiara Brust

Date: 8/31/19

Goal:

1. Glycerol stocks of mCherry from 8/29/19 transformation

Protocol:

Glycerol Stocks

1. Took 1 mL of 50% glycerol and 1 ml of the overnight culture (after incubation) and added to a glycerol stock tube.
2. Labeled with name, date, and the contents and store in the -80° C freezer in CLSO 442